

Maxer  
GERMANY

ARTHROSCOPY



**M** ARTHROSCOPY  
**SYSTEM**

[maxerendoscopy.com](http://maxerendoscopy.com)



## ISO Certificate

Quality management system as per  
ISO 13485 :2012 Medical Device Certification  
for development & Production & Endoscopic  
instrument & Equipments

## CE Certificate

As per council directive 93/42/EEC  
Concerning medical devices



## About Maxer

Maxer Medizintechnik GmbH has established itself as a symbol of creativity and latest technology and is an internationally qualified supplier of world class medical instruments and devices.

Specialized in the field of Minimal Invasive Surgery, Maxer manufactures wide range of equipments such as Laparoscopy, Hysteroscopy, Urology, Pediatric Surgery, Arthroscopy, Spine Endoscopy, ENT etc.

Maxer has been one of the foremost companies to bring out HD technology in MIS field. In 2004, Maxer brought out first endoscopy camera with Matrix Metering. In 2008, Maxer brought out HD Telescopes. In 2009, Maxer launched first Networked Endoscopy Camera System. In 2010, Maxer introduced pathbreaking 16 Fr Resectoscope for Office Hysteroscopy. In 2011, Maxer launched LED Headlight with 200,000 Lux illumination, again a state of art in this segment. In 2012, Maxer launched LED Light sources to keep Maxer customers always at the cutting edge of technology. Maxer has also launched HD cameras with integrated Full HD recording. 3 CMOS HD cameras with 6 megapixel camera head will be launched soon.

Owners of Maxer have been working in Medical market for more than 20 years and have installed thousands of sets with high level of customer satisfaction. This background has given Maxer profound credibility in the field of endoscopy.

Due to vibrant co-operation with end users, Maxer is able to continuously improve upon the product range and also able to meet special modification demands. Maxer believes that growth and innovation go hand-in-hand.









Maxer sells complete system with a focus on product quality and customers satisfaction. Maxer provides sales and service training to it's distributor in order to be successful.

Maxer believes that sound ethics with professionalism will only give success to the company. Maxer's growth and success are based on long term partnerships with customers, distributors. Maxer currently exports it's products to Europe, Asia, Africa.

Maxer is headed by Mr. Giovanni Miranda in the capacity as a Managing Director.

How to reach to Maxer - Maxer is located in Würlingen (sales office) within the medical instruments manufacturing region around Tuttlingen in South Germany. Maxer also has registered office at Spaichingen.




## INDEX

|                                                 |            |
|-------------------------------------------------|------------|
| Arthroscope & Information                       | A-4 to 5   |
| Arthroscopy Sheaths                             | A-6 to 7   |
| Tube Shaft Instruments                          | A-8 to 15  |
| Max-Line (Pinless Hinge) Instrument Information | A-16 to 18 |
| Max-Line Knee Instruments                       | A-19 to 34 |
| Max-Line Small Joint Instruments                | A-35 to 37 |
| Max-Line Handles for Shoulder Instruments       | A-38       |
| Max-Line Shoulder Instruments                   | A-39 to 42 |
| Max-Line Hip Instruments                        | A-43 to 46 |
| Arthroscopy Probes                              | A-47 to 48 |
| Shaver Blades                                   | A-49 to 51 |
| Monopolar Arthroscopy Electrodes                | A-52 to 54 |
| Bipolar Arthroscopy Electrodes                  | A-55 to 56 |

## Arthroscopes from Maxer, Germany

Maxer has been setting innovative standards for many years. Maxer laser welded Arthroscope can withstand numerous autoclaving cycles. Its glass cone technology is also unique. Maxer endoscopes are known for third generation optics which gives high resolution and brightness all over the image.



### Sapphire tip -

For 3 times more scratch resistance,  
safe with shaver blades

Ø 4.0 mm




Ø 2.7 mm



### Laser welded joints -


No chance of leakage  
even after numerous  
autoclaving cycle



### Glass cone -

For longer life of fibers

Complete stainless steel -  
Increased sturdiness



## FEATURES

- High Definition Laparoscopes have better optics to give enhanced clarity & brightness
- High-quality, multi coated (anti-reflective) optical system
- Highly refractive glass for optimal color rendering even when the image is magnified
- Extremely precise centering of the optical and mechanical components for an image that is brilliant and high contrast even around the edges
- Triple tube design for higher stability with many types of endoscopes
- Nitrogen filled optical system to avoid interior condensation
- High resolution scopes, Sapphire cover glass on both ends,
- Glass cone to collect light from F. O. Cable, Autoclavable

## Arthroscopes



### Standard Version 4 mm

| 0°                | 30°               | 45°               | 70°               | Working Length | Ø mm |
|-------------------|-------------------|-------------------|-------------------|----------------|------|
| <b>15.20.4000</b> | <b>15.20.4030</b> | <b>15.20.4045</b> | <b>15.20.4070</b> | 185 mm         | 4.0  |

### HD Version 4 mm

| 0°                | 30°               | 45°               | 70°               | Working Length | Ø mm |
|-------------------|-------------------|-------------------|-------------------|----------------|------|
| <b>15.20.4100</b> | <b>15.20.4130</b> | <b>15.20.4145</b> | <b>15.20.4170</b> | 185 mm         | 4.0  |

### Standard Version 2.7 mm


| 0°                | 30°               | 45° | 70°               | Working Length | Ø mm |
|-------------------|-------------------|-----|-------------------|----------------|------|
| <b>15.20.2700</b> | <b>15.20.2730</b> | -   | <b>15.20.2770</b> | 195 mm         | 2.7  |
| <b>15.20.2900</b> | <b>15.20.2930</b> | -   | <b>15.20.2970</b> | 140 mm         | 2.7  |

### Standard Version 1.9 mm

| 0°                | 30°               | 45° | 70° | Working Length | Ø mm |
|-------------------|-------------------|-----|-----|----------------|------|
| <b>15.20.1900</b> | <b>15.20.1930</b> | -   | -   | 60 mm          | 1.9  |
| <b>15.20.1910</b> | <b>15.20.1940</b> | -   | -   | 100 mm         | 1.9  |

## Arthroscopy Sheath

Easy Locking System



For Arthroscopes 4 mm, working length 185 mm

|                     |                                                         |
|---------------------|---------------------------------------------------------|
| <b>50.80.4030EL</b> | High flow trocar sheath 6 mm with 2 rotating stopcocks, |
| <b>50.80.4031</b>   | Obturator blunt 5 mm                                    |
| <b>50.80.4032</b>   | Obturator conical blunt 5 mm                            |
| <b>50.80.4033</b>   | Obturator pyramidal (sharp) 5 mm                        |



## Arthroscopy Sheath

For Arthroscopes 2.7 mm, working length 195 mm

|                   |                                     |
|-------------------|-------------------------------------|
| <b>50.80.2730</b> | Trocar sheath with 1 fixed stopcock |
| <b>50.80.2731</b> | Trocar blunt                        |
| <b>50.80.2732</b> | Trocar conical blunt                |
| <b>50.80.2733</b> | Trocar pyramidal (sharp)            |

For Arthroscopes 2.7 mm, working length 140 mm

|                   |                                     |
|-------------------|-------------------------------------|
| <b>50.80.2930</b> | Trocar sheath with 1 fixed stopcock |
| <b>50.80.2934</b> | Trocar sheath with 2 fixed stopcock |
| <b>50.80.2931</b> | Obturator blunt                     |
| <b>50.80.2932</b> | Obturator conical blunt             |
| <b>50.80.2933</b> | Obturator pyramidal (sharp)         |

For Arthroscopes 1.9 mm, working length 60 mm

|                   |                                     |
|-------------------|-------------------------------------|
| <b>50.80.1930</b> | Trocar sheath with 1 fixed stopcock |
| <b>50.80.1935</b> | Trocar sheath with 2 fixed stopcock |
| <b>50.80.1931</b> | Obturator blunt                     |
| <b>50.80.1932</b> | Obturator conical blunt             |
| <b>50.80.1933</b> | Obturator pyramidal (sharp)         |

## Arthroscopy Instruments

Standard Instruments



### FEATURES

- Tubular Construction
- Built in overload Protection
- Stainless Steel
- Autoclavable 134° C

|   |                 |
|---|-----------------|
| ■ | straight        |
| ▲ | curved upward   |
| ▼ | curved downward |
| ◀ | curved left     |
| ▶ | curved right    |

## Arthroscopy Instruments



Hook punch, curved vertical, without spoon

|              | Ø 3.4 mm   | Ø 2.8 mm   |   |     |
|--------------|------------|------------|---|-----|
| 15° upward   | 50.50.1534 | 50.50.1528 | ▲ | 15° |
| 7° upward    | 50.50.0734 | 50.50.0728 | ▲ | 7°  |
| straight     | 50.50.3400 | 50.50.2800 | ■ | 0°  |
| 7° downward  | 50.50.3407 | 50.50.2807 | ▼ | 7°  |
| 15° downward | 50.50.3415 | 50.50.2815 | ▼ | 15° |

Curved horizontal

|           | Ø 3.4 mm   | Ø 2.8 mm   |   |     |
|-----------|------------|------------|---|-----|
| 45° right | 50.50.4534 | 50.50.4528 | ► | 45° |
| 30° right | 50.50.3034 | 50.50.3028 | ► | 30° |
| 30° left  | 50.50.3430 | 50.50.2830 | ◄ | 30° |
| 45° left  | 50.50.3445 | 50.50.2845 | ◄ | 45° |

■ straight ▲ curved upward ▼ curved downward ◄ curved left ► curved right

## Arthroscopy Instruments



Hook scissor, Curved Vertical

|              | Ø 3.4 mm   | Ø 2.8 mm   |   |     |
|--------------|------------|------------|---|-----|
| 15° upward   | 50.34.1500 | 50.28.1500 | ▲ | 15° |
| 7° upward    | 50.34.0700 | 50.28.0700 | ▲ | 7°  |
| straight     | 50.34.0000 | 50.28.0000 | ■ | 0°  |
| 7° downward  | 50.34.0007 | 50.28.0007 | ▼ | 7°  |
| 15° downward | 50.34.0015 | 50.28.0015 | ▼ | 15° |

Curved horizontal

|           | Ø 3.4 mm   | Ø 2.8 mm   |   |     |
|-----------|------------|------------|---|-----|
| 45° right | 50.34.4500 | 50.28.4500 | ► | 45° |
| 30° right | 50.34.3000 | 50.28.3000 | ► | 30° |
| 30° left  | 50.34.0030 | 50.28.0030 | ◄ | 30° |
| 45° left  | 50.34.0045 | 50.28.0045 | ◄ | 45° |

■ straight ▲ curved upward ▼ curved downward ◄ curved left ► curved right



## Arthroscopy Instruments



Punch slim, curved vertical

|              | Ø 3.4 mm   | Ø 2.8 mm   |   |     |
|--------------|------------|------------|---|-----|
| 15° upward   | 50.60.1534 | 50.60.1528 | ▲ | 15° |
| 7° upward    | 50.60.0734 | 50.60.0728 | ▲ | 7°  |
| straight     | 50.60.3400 | 50.60.2800 | ■ | 0°  |
| 7° downward  | 50.60.3407 | 50.60.2807 | ▼ | 7°  |
| 15° downward | 50.60.3415 | 50.60.2815 | ▼ | 15° |

Curved horizontal

|           | Ø 3.4 mm   | Ø 2.8 mm   |   |     |
|-----------|------------|------------|---|-----|
| 45° right | 50.60.4534 | 50.60.4528 | ► | 45° |
| 30° right | 50.60.3034 | 50.60.3028 | ► | 30° |
| 30° left  | 50.60.3430 | 50.60.2830 | ◄ | 30° |
| 45° left  | 50.60.3445 | 50.60.2845 | ◄ | 45° |

■ straight ▲ curved upward ▼ curved downward ◄ curved left ► curved right

## Arthroscopy Instruments



Square punch wide, curved vertical

|  |            |   |     |
|--|------------|---|-----|
|  | Ø 3.4 mm   |   |     |
|  | 50.70.1534 | ▲ | 15° |
|  | 50.70.0734 | ▲ | 7°  |
|  | 50.70.3400 | ■ | 0°  |
|  | 50.70.3407 | ▼ | 7°  |
|  | 50.70.3415 | ▼ | 15° |

Curved horizontal


|  |            |   |     |
|--|------------|---|-----|
|  | Ø 3.4 mm   |   |     |
|  | 50.70.4534 | ► | 45° |
|  | 50.70.3034 | ► | 30° |
|  | 50.70.3430 | ◄ | 30° |
|  | 50.70.3445 | ◄ | 45° |

■ straight ▲ curved upward ▼ curved downward ◄ curved left ► curved right

## Arthroscopy Instruments




Square punch wide, 90°, curved vertical, Ø 3.4 mm



| 90° right         | 90° left          |   |     |
|-------------------|-------------------|---|-----|
| <b>50.70.1590</b> | <b>50.75.1590</b> | ▲ | 15° |
| <b>50.70.0790</b> | <b>50.75.0790</b> | ▲ | 7°  |
| <b>50.70.9000</b> | <b>50.75.9000</b> | ■ | 0°  |
| <b>50.70.9007</b> | <b>50.75.9007</b> | ▼ | 7°  |
| <b>50.70.9015</b> | <b>50.75.9015</b> | ▼ | 15° |

Curved horizontal



| 90° right         | 90° left          |   |     |
|-------------------|-------------------|---|-----|
| <b>50.70.4590</b> | <b>50.75.4590</b> | ► | 45° |
| <b>50.70.3090</b> | <b>50.75.3090</b> | ► | 30° |
| <b>50.70.9030</b> | <b>50.75.9030</b> | ◄ | 30° |
| <b>50.70.9045</b> | <b>50.75.9045</b> | ◄ | 45° |

■ straight ▲ curved upward ▼ curved downward ◄ curved left ► curved right



## Arthroscopy Instruments



Probe excisions & reduction forceps, with thorn, straight tube

| without spreading spring | with spreading spring |          |
|--------------------------|-----------------------|----------|
| <b>50.45.4500</b>        | <b>50.46.4500</b>     | Ø 4.5 mm |



“Alligator” grasping forceps, straight tube

| without ratchet   | with ratchet      |          |
|-------------------|-------------------|----------|
| <b>50.45.4528</b> | <b>50.46.4528</b> | Ø 2.8 mm |
| <b>50.45.4534</b> | <b>50.46.4534</b> | Ø 3.4 mm |
| <b>50.45.4545</b> | <b>50.46.4545</b> | Ø 4.5 mm |



Grasping forceps, with 1 x 2 teeth, straight tube

| without ratchet   | with ratchet      |          |
|-------------------|-------------------|----------|
| <b>50.45.3401</b> | <b>50.46.3401</b> | Ø 3.4 mm |



Biopsy & grasping forceps, straight tube


| without ratchet   | with ratchet      |          |
|-------------------|-------------------|----------|
| <b>50.45.3402</b> | <b>50.46.3402</b> | Ø 3.4 mm |



|                   |                                     |
|-------------------|-------------------------------------|
| <b>50.46.3405</b> | Straight shaft grasper with ratchet |
|-------------------|-------------------------------------|

## Max-Line Arthroscopy Instruments

Pinless Hinge



- Newly designed using pin - less hinge that allows application of higher force and superior cutting efficiency.
- New design also facilitates smaller tip profile enabling surgeon to operate in harder to reach areas.
- Slotted single piece outer shaft design provides sturdy design and longer durability.
- A wide range of tips is available to cover broad spectrum of treatments.
- All the instruments are made in Germany using finest craftsmanship.
- Most Max-Line instruments are Pinless.

## Max-Line Arthroscopy Instruments

Knee Arthroscopy



Tip Profile



Bite Width



Tip Width



Straight

The primary style for direct approach



Angled 1/4 curve at tip increased access to the posterior meniscus




Angled shaft follows the contour of the condyle and keeps the cutting surface parallel to the tibial plateau



Angled to reach the anterior and posterior meniscus





## Max-Line Arthroscopy Instruments

Handles



Standards handles

**50.80.----**

Scissors, basket punches and graspers feature our traditional Standard handles



Standard handle with ratchet

**50.83.----**

All graspers featured this infinite-position, locking-handle design

Cigar handles

**50.81.----**

Available on rotary basket punches and scissors



Roratable handles

**50.82.----**

Available on rotary basket punches and scissors

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles

Hook Punch, slim style



| Standard Handle   | Description           | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------------|-------------|------------|-----------|
| <b>50.80.5100</b> | upbiter, straight     | 1.90 mm     | 1.67 mm    | 2.89 mm   |
| <b>50.80.5101</b> | upbiter, curved left  | 1.90 mm     | 1.67 mm    | 2.89 mm   |
| <b>50.80.5102</b> | upbiter, curved right | 1.90 mm     | 1.67 mm    | 2.89 mm   |

Hook Punch, blunt nose 1.0mm



| Standard Handle   | Description  | Tip Profile | Bite Width | Tip Width |
|-------------------|--------------|-------------|------------|-----------|
| <b>50.80.5103</b> | straight     | 1.83 mm     | 1.67 mm    | 3.32 mm   |
| <b>50.80.5104</b> | upswept      | 1.83 mm     | 1.67 mm    | 3.32 mm   |
| <b>50.80.5105</b> | curved left  | 1.83 mm     | 1.67 mm    | 3.32 mm   |
| <b>50.80.5106</b> | curved right | 1.83 mm     | 1.67 mm    | 3.32 mm   |

Hook Punch, short cut



| Standard Handle   | Description |
|-------------------|-------------|
| <b>50.80.5107</b> | straight    |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Duckbill Basket Punch, small

| Standard Handle   | Description | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------|-------------|------------|-----------|
| <b>50.80.5111</b> | straight    | 2.52 mm     | 2.50 mm    | 4.15 mm   |



Duckbill Basket Punch, small

| Standard Handle   | Description       | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------------|-------------|------------|-----------|
| <b>50.80.5112</b> | upbiter, straight | 2.52 mm     | 2.50 mm    | 4.15 mm   |



Duckbill Basket Punch, medium

| Standard Handle   | Description  | Tip Profile | Bite Width | Tip Width |
|-------------------|--------------|-------------|------------|-----------|
| <b>50.80.5121</b> | straight     | 2.52 mm     | 3.20 mm    | 5.05 mm   |
| <b>50.80.5122</b> | curved left  | 2.52 mm     | 3.20 mm    | 5.05 mm   |
| <b>50.80.5123</b> | curved right | 2.52 mm     | 3.20 mm    | 5.05 mm   |
| <b>50.80.5124</b> | upswept      | 2.52 mm     | 3.20 mm    | 5.05 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Duckbill Basket Punch, medium

| Standard Handle | Description           | Tip Profile | Bite Width | Tip Width |
|-----------------|-----------------------|-------------|------------|-----------|
| 50.80.5125      | upbiter, straight     | 2.52 mm     | 3.20 mm    | 5.05 mm   |
| 50.80.5126      | upbiter, curved left  | 2.52 mm     | 3.20 mm    | 5.05 mm   |
| 50.80.5127      | upbiter, curved right | 2.52 mm     | 3.20 mm    | 5.05 mm   |
| 50.80.5128      | upbiter, upswept      | 2.52 mm     | 3.20 mm    | 5.05 mm   |



Duckbill Basket Punch, large

| Standard Handle | Description | Tip Profile | Bite Width | Tip Width |
|-----------------|-------------|-------------|------------|-----------|
| 50.80.5131      | straight    | 2.52 mm     | 4.50 mm    | 6.15 mm   |



Duckbill Basket Punch, large

| Standard Handle | Description       | Tip Profile | Bite Width | Tip Width |
|-----------------|-------------------|-------------|------------|-----------|
| 50.80.5132      | upbiter, straight | 2.52 mm     | 4.50 mm    | 6.15 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Duckbill Basket Punch

| Standard Handle   | Description                | Tip Profile | Bite Width | Tip Width |
|-------------------|----------------------------|-------------|------------|-----------|
| <b>50.80.5136</b> | diagonal cut, curved left  | 2.52 mm     | 3.20 mm    | 5.05 mm   |
| <b>50.80.5137</b> | diagonal cut, curved right | 2.52 mm     | 3.20 mm    | 5.05 mm   |



Duckly Basket Punch, small

| Standard Handle   | Description  | Tip Profile | Bite Width | Tip Width |
|-------------------|--------------|-------------|------------|-----------|
| <b>50.80.5144</b> | straight     | 1.93 mm     | 2.50 mm    | 4.15 mm   |
| <b>50.80.5145</b> | curved left  | 1.93 mm     | 2.50 mm    | 4.15 mm   |
| <b>50.80.5146</b> | curved right | 1.93 mm     | 2.50 mm    | 4.15 mm   |
| <b>50.80.5147</b> | upswept      | 1.93 mm     | 2.50 mm    | 4.15 mm   |



Duckly Basket Punch, small

| Standard Handle   | Description       | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------------|-------------|------------|-----------|
| <b>50.80.5148</b> | Upbiter, straight | 1.93 mm     | 2.50 mm    | 4.15 mm   |



## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Duckly Basket Punch, medium

| Standard Handle   | Description  | Tip Profile | Bite Width | Tip Width |
|-------------------|--------------|-------------|------------|-----------|
| <b>50.80.5155</b> | straight     | 1.93 mm     | 3.20 mm    | 5.00 mm   |
| <b>50.80.5156</b> | curved left  | 1.93 mm     | 3.20 mm    | 5.00 mm   |
| <b>50.80.5157</b> | curved right | 1.93 mm     | 3.20 mm    | 5.00 mm   |
| <b>50.80.5158</b> | upswept      | 1.93 mm     | 3.20 mm    | 5.00 mm   |

Duckly Basket Punch, medium

| Standard Handle   | Description           | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------------|-------------|------------|-----------|
| <b>50.80.5159</b> | upbiter, straight     | 1.93 mm     | 3.20 mm    | 5.00 mm   |
| <b>50.80.5160</b> | upbiter, curved left  | 1.93 mm     | 3.20 mm    | 5.00 mm   |
| <b>50.80.5161</b> | upbiter, curved right | 1.93 mm     | 3.20 mm    | 5.00 mm   |
| <b>50.80.5162</b> | upbiter, upswept      | 1.93 mm     | 3.20 mm    | 5.00 mm   |

Duckly Basket Punch, large

| Standard Handle   | Description | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------|-------------|------------|-----------|
| <b>50.80.5166</b> | straight    | 1.93 mm     | 4.50 mm    | 6.15 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Duckly Basket Punch, large

| Standard Handle | Description       | Tip Profile | Bite Width | Tip Width |
|-----------------|-------------------|-------------|------------|-----------|
| 50.80.5167      | upbiter, straight | 1.93 mm     | 4.50 mm    | 6.15 mm   |



Duckly Punch Elevator

| Standard Handle | Description       | Tip Profile | Bite Width | Tip Width |
|-----------------|-------------------|-------------|------------|-----------|
| 50.80.5171      | upbiter, straight | 1.93 mm     | 3.20 mm    | 5.00 mm   |



Basket Punch, slim style

| Standard Handle | Description | Tip Profile | Bite Width | Tip Width |
|-----------------|-------------|-------------|------------|-----------|
| 50.80.5176      | straight    | 1.90 mm     | 1.67 mm    | 2.89 mm   |
| 50.80.5177      | upswept     | 1.90 mm     | 1.67 mm    | 2.89 mm   |



Basket Punch, slim style

| Standard Handle | Description          | Tip Profile | Bite Width | Tip Width |
|-----------------|----------------------|-------------|------------|-----------|
| 50.80.5178      | jaw 45° angled left  | 1.50 mm     | 1.67 mm    | 2.89 mm   |
| 50.80.5179      | jaw 45° angled right | 1.50 mm     | 1.67 mm    | 2.89 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Basket Punch, oval

| Standard Handle   | Description  | Tip Profile | Bite Width | Tip Width |
|-------------------|--------------|-------------|------------|-----------|
| <b>50.80.5185</b> | straight     | 1.65 mm     | 4.39 mm    | 6.57 mm   |
| <b>50.80.5186</b> | curved left  | 1.65 mm     | 4.39 mm    | 6.57 mm   |
| <b>50.80.5187</b> | curved right | 1.65 mm     | 4.39 mm    | 6.57 mm   |
| <b>50.80.5188</b> | upswept      | 1.65 mm     | 4.39 mm    | 6.57 mm   |

Basket Punch, oval

| Standard Handle   | Description           | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------------|-------------|------------|-----------|
| <b>50.80.5189</b> | upbiter, straight     | 1.65 mm     | 4.39 mm    | 6.57 mm   |
| <b>50.80.5190</b> | upbiter, curved left  | 1.65 mm     | 4.39 mm    | 6.57 mm   |
| <b>50.80.5191</b> | upbiter, curved right | 1.65 mm     | 4.39 mm    | 6.57 mm   |
| <b>50.80.5192</b> | upbiter, upswept      | 1.65 mm     | 4.39 mm    | 6.57 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Basket Punch spoon

| Standard Handle   | Description          | Tip Profile | Bite Width | Tip Width |
|-------------------|----------------------|-------------|------------|-----------|
| <b>50.80.5195</b> | 1.10 mm, straight    | 1.10 mm     | 1.59 mm    | 3.25 mm   |
| <b>50.80.5196</b> | 1.10 mm, upswept     | 1.10 mm     | 1.59 mm    | 3.25 mm   |
| <b>50.80.5201</b> | 1.3 mm, curved left  | 1.30 mm     | 1.59 mm    | 3.25 mm   |
| <b>50.80.5202</b> | 1.3 mm, curved right | 1.30 mm     | 1.59 mm    | 3.25 mm   |

Basket Punch spoon, 1.3 mm

| Standard Handle   | Description           | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------------|-------------|------------|-----------|
| <b>50.80.5203</b> | upbiter, straight     | 1.60 mm     | 1.59 mm    | 3.27 mm   |
| <b>50.80.5204</b> | upbiter, curved left  | 1.60 mm     | 1.59 mm    | 3.27 mm   |
| <b>50.80.5205</b> | upbiter, curved right | 1.60 mm     | 1.59 mm    | 3.27 mm   |

Basket Punch spoon

| Standard Handle   | Description      | Tip Profile | Bite Width | Tip Width |
|-------------------|------------------|-------------|------------|-----------|
| <b>50.80.5206</b> | 1.3 mm, upswept  | 1.30 mm     | 1.59 mm    | 3.27 mm   |
| <b>50.80.5207</b> | 1.5 mm, straight | 1.60 mm     | 1.59 mm    | 3.88 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles

Basket Punch spoon



| Standard Handle   | Description       | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------------|-------------|------------|-----------|
| <b>50.80.5208</b> | 1.3 mm, downbiter | 1.30 mm     | 1.59 mm    | 3.27 mm   |

Basket Punch spoon



| Standard Handle   | Description      | Tip Profile | Bite Width | Tip Width |
|-------------------|------------------|-------------|------------|-----------|
| <b>50.80.5209</b> | 1.5 mm, upswept  | 1.60 mm     | 1.59 mm    | 3.88 mm   |
| <b>50.80.5210</b> | 2.7 mm, straight | 2.50 mm     | 1.77 mm    | 3.78 mm   |
| <b>50.80.5211</b> | 3.4 mm, straight | 3.00 mm     | 2.03 mm    | 3.76 mm   |

Basket Punch spoon



| Standard Handle   | Description       | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------------|-------------|------------|-----------|
| <b>50.80.5212</b> | upbiter, straight | 2.28 mm     | 1.59 mm    | 3.88 mm   |

Basket Punch meniscal elevator



| Standard Handle   | Description      | Tip Profile | Bite Width | Tip Width |
|-------------------|------------------|-------------|------------|-----------|
| <b>50.80.5216</b> | 1.3 mm, straight | 1.90 mm     | 1.60 mm    | 3.20 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Basket Punch meniscal elevator

| Standard Handle   | Description               | Tip Profile | Bite Width | Tip Width |
|-------------------|---------------------------|-------------|------------|-----------|
| <b>50.80.5217</b> | 1.3 mm, upbiter, straight | 1.90 mm     | 1.60 mm    | 3.20 mm   |



Posterior Basket Punch

| Standard Handle   | Description  | Tip Profile | Bite Width | Tip Width |
|-------------------|--------------|-------------|------------|-----------|
| <b>50.80.5221</b> | straight     | 2.46 mm     | 2.18 mm    | 4.00 mm   |
| <b>50.80.5222</b> | curved left  | 2.46 mm     | 2.18 mm    | 4.00 mm   |
| <b>50.80.5223</b> | curved right | 2.46 mm     | 2.18 mm    | 4.00 mm   |
| <b>50.80.5224</b> | upswept      | 2.46 mm     | 2.18 mm    | 4.00 mm   |



Posterior Basket Punch

| Standard Handle   | Description           | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------------|-------------|------------|-----------|
| <b>50.80.5225</b> | upbiter, straight     | 2.46 mm     | 2.18 mm    | 4.00 mm   |
| <b>50.80.5226</b> | upbiter, curved left  | 2.46 mm     | 2.18 mm    | 4.00 mm   |
| <b>50.80.5227</b> | upbiter, curved right | 2.46 mm     | 2.18 mm    | 4.00 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Basket Punch unicorn

| Standard Handle   | Description     | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------|-------------|------------|-----------|
| <b>50.80.5231</b> | straight, left  | 3.93 mm     | 2.38 mm    | 5.58 mm   |
| <b>50.80.5232</b> | straight, right | 3.93 mm     | 2.38 mm    | 5.58 mm   |



Rongeur Raptor

| Standard Handle   | Description | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------|-------------|------------|-----------|
| <b>50.80.5233</b> | straight    | -           | -          | -         |



Scissors Punch

| Standard Handle   | Description  | Tip Profile | Bite Width | Tip Width |
|-------------------|--------------|-------------|------------|-----------|
| <b>50.80.5236</b> | straight     | 2.20 mm     | 1.20 mm    | 2.80 mm   |
| <b>50.80.5237</b> | curved left  | 2.20 mm     | 1.20 mm    | 2.80 mm   |
| <b>50.80.5238</b> | curved right | 2.20 mm     | 1.20 mm    | 2.80 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Standards handles



Scissors Punch

| Standard Handle   | Description           | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------------|-------------|------------|-----------|
| <b>50.80.5236</b> | upbiter, straight     | 2.20 mm     | 1.20 mm    | 2.80 mm   |
| <b>50.80.5237</b> | upbiter, curved left  | 2.20 mm     | 1.20 mm    | 2.80 mm   |
| <b>50.80.5238</b> | upbiter, curved right | 2.20 mm     | 1.20 mm    | 2.80 mm   |

Hook Scissors

| Standard Handle   | Description           | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------------|-------------|------------|-----------|
| <b>50.80.5239</b> | upbiter, straight     | 2.20 mm     | 1.20 mm    | 2.80 mm   |
| <b>50.80.5240</b> | upbiter, curved left  | 2.20 mm     | 1.20 mm    | 2.80 mm   |
| <b>50.80.5241</b> | upbiter, curved right | 2.20 mm     | 1.20 mm    | 2.80 mm   |



## Max-Line Arthroscopy Instruments

Knee Instruments, Punches



Rotatable handles



Cigar Handle



### Rotary Scissors

| Rotatable Handle | Cigar Handle | Description           |
|------------------|--------------|-----------------------|
| 50.82.5251       | 50.81.5251   | jaw 20°, angled left  |
| 50.82.5252       | 50.81.5252   | jaw 20°, angled right |

### Rotary Scissors, serrated

| Rotatable Handle | Cigar Handle | Description           |
|------------------|--------------|-----------------------|
| 50.82.5253       | 50.81.5253   | jaw 20°, angled left  |
| 50.82.5254       | 50.81.5254   | jaw 20°, angled right |

### Rotary Hook Scissors

| Rotatable Handle | Cigar Handle | Description           |
|------------------|--------------|-----------------------|
| 50.82.5255       | 50.81.5255   | jaw 60°, angled left  |
| 50.82.5256       | 50.81.5256   | jaw 60°, angled right |

## Max-Line Arthroscopy Instruments

Knee Instruments, Rotary Basket Punches, Grasping Forceps



Rotatable handles



Cigar Handle



### Rotary Basket Punch

| Rotatable Handle  | Cigar Handle      | Description                   | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------------|-------------------------------|-------------|------------|-----------|
| <b>50.82.5300</b> | <b>50.81.5300</b> | 2.2 mm, jaw 90°, angled left  | 1.60 mm     | 2.20 mm    | 4.10 mm   |
| <b>50.82.5301</b> | <b>50.81.5301</b> | 2.2 mm, jaw 90°, angled right | 1.60 mm     | 2.20 mm    | 4.10 mm   |
| <b>50.82.5302</b> | <b>50.81.5302</b> | 3.4 mm, jaw 90°, angled left  | 1.60 mm     | 3.40 mm    | 5.30 mm   |
| <b>50.82.5303</b> | <b>50.81.5303</b> | 3.4 mm, jaw 90°, angled right | 1.60 mm     | 3.40 mm    | 5.30 mm   |

## Max-Line Arthroscopy Instruments

Knee Instruments, Rotary Basket Punches, Grasping Forceps



Standard handle with ratchet

### Grasping Forceps

| Standard handle with ratchet | Description         |
|------------------------------|---------------------|
| <b>50.83.5400</b>            | 1x2 teeth, straight |



### Grasping Forceps cupped, serrated

| Standard handle with ratchet | Description      |
|------------------------------|------------------|
| <b>50.83.5401</b>            | 2.7 mm, straight |
| <b>50.83.5402</b>            | 3.4 mm, straight |



### Grasping Forceps Alligator, 2.7 mm

| Standard handle with ratchet | Description       |
|------------------------------|-------------------|
| <b>50.83.5403</b>            | upbiter, straight |
| <b>50.83.5404</b>            | upbiter, upswept  |

## Max-Line Arthroscopy Instruments

Knee Instruments, Grasping Forceps



Standard handle with ratchet

Grasping Forceps Alligator, 2.7 mm



| Standard handle with ratchet | Description |
|------------------------------|-------------|
| 50.83.5405                   | straight    |
| 50.83.5406                   | upswept     |

Grasping Forceps Alligator, 3.4 mm



| Standard handle with ratchet | Description |
|------------------------------|-------------|
| 50.83.5407                   | straight    |
| 50.83.5408                   | upswept     |



Tissue Grasping Forceps with channel 2 mm



| Standard handle with ratchet | Description |
|------------------------------|-------------|
| 50.83.5416                   | straight    |

## Max-Line Arthroscopy Instruments

Small Joint-Instruments,Punches



Standards handles



Punch, Oval Jr.

| Standard Handle   | Description  | Tip Profile | Bite Width | Tip Width |
|-------------------|--------------|-------------|------------|-----------|
| <b>50.80.5266</b> | straight     | 1.20 mm     | 2.30 mm    | 3.50 mm   |
| <b>50.80.5267</b> | curved left  | 1.20 mm     | 2.30 mm    | 3.50 mm   |
| <b>50.80.5268</b> | curved right | 1.20 mm     | 2.30 mm    | 3.50 mm   |



Duckbill-Punch Jr.

| Standard Handle   | Description | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------|-------------|------------|-----------|
| <b>50.80.5271</b> | straight    | 1.20 mm     | 2.30 mm    | 3.50 mm   |




Duckbill-Punch Jr.

| Standard Handle   | Description       | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------------|-------------|------------|-----------|
| <b>50.80.5272</b> | upbiter, straight | 1.20 mm     | 2.30 mm    | 3.50 mm   |



Basket-Punch Blunt, Jr

| Standard Handle   | Description | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------|-------------|------------|-----------|
| <b>50.80.5275</b> | straight    | 1.70 mm     | 1.80 mm    | 3.00 mm   |



**Max-Line Arthroscopy Instruments**

Small Joint-Instruments, Scissors, Rongeurs, Grasping Forceps



Standards handles

Hook Scissors Jr.



| Standard Handle | Description |
|-----------------|-------------|
| 50.80.5278      | straight    |

Rongeur Raptor Jr.



| Standard Handle | Description |
|-----------------|-------------|
| 50.80.5280      | straight    |

## Max-Line Arthroscopy Instruments

Small Joint-Instruments, Scissors, Rongeurs, Grasping Forceps



Standards handles

Standard handle with ratchet



Grasping Forceps Alligator Jr.

| Standard handle | Standard handle with ratchet | Description |
|-----------------|------------------------------|-------------|
| -               | 50.83.5283                   | straight    |



Grasping Forceps Alligator Jr., 1x2 teeth

| Standard handle | Standard handle with ratchet | Description |
|-----------------|------------------------------|-------------|
| -               | 50.83.5284                   | straight    |

## Max-Line Arthroscopy Instruments

### Handles



#### Standard Handle

**50.80.----**

Scissors, basket punches and graspers feature our traditional standard handles



#### Standard handle with ratchet

**50.83.----**

All graspers featured this infinite-position, locking-handle design



#### Axial Handle

**50.85.----**

All graspers featured this infinite-position, handle design



## Max-Line Arthroscopy Instruments

Shoulder-Instruments, Suture Guide, Suture Hook



Standard Handle



Axial Handle



Clever Hook, Suture Hook

| Standard Handle   | Axial Handle      | Description                 | Working Length |
|-------------------|-------------------|-----------------------------|----------------|
| <b>50.80.5500</b> | <b>50.85.5500</b> | pointed tip left, straight  | 155 mm         |
| <b>50.80.5501</b> | <b>50.85.5501</b> | pointed tip right, straight | 155 mm         |



Sixer, Suture Hook

| Standard Handle   | Axial Handle      | Description  | Working Length |
|-------------------|-------------------|--------------|----------------|
| <b>50.80.5504</b> | <b>50.85.5504</b> | turned left  | 155 mm         |
| <b>50.80.5505</b> | <b>50.85.5505</b> | turned right | 155 mm         |



Bird Beak, Suture Hook Forceps

| Standard Handle   | Axial Handle      | Description | Working Length |
|-------------------|-------------------|-------------|----------------|
| <b>50.80.5508</b> | <b>50.85.5508</b> | straight    | 155 mm         |



Bird Beak, Suture Hook Forceps

| Standard Handle   | Axial Handle      | Description                    | Working Length |
|-------------------|-------------------|--------------------------------|----------------|
| <b>50.80.5509</b> | <b>50.85.5509</b> | jaw 45° angled left, straight  | 155 mm         |
| <b>50.80.5510</b> | <b>50.85.5510</b> | jaw 45° angled right, straight | 155 mm         |

## Max-Line Arthroscopy Instruments

Shoulder-Instruments, Suture Guide, Suture Hook



Standard Handle

Axial Handle



Bird Beak, Suture Hook Forceps

| Standard Handle   | Axial Handle      | Description                       | Working Length |
|-------------------|-------------------|-----------------------------------|----------------|
| <b>50.80.5511</b> | <b>50.85.5511</b> | jaw 20° angled, upbiter, straight | 155 mm         |
| <b>50.80.5512</b> | <b>50.85.5512</b> | jaw 40° angled, upbiter, straight | 155 mm         |



Suture Punch

| Standard Handle   | Axial Handle      | Description | Working Length |
|-------------------|-------------------|-------------|----------------|
| <b>50.80.5625</b> | <b>50.85.5625</b> | straight    | 155 mm         |



Double Hook, Suture Hook

| Standard Handle   | Axial Handle      | Description | Working Length |
|-------------------|-------------------|-------------|----------------|
| <b>50.80.5626</b> | <b>50.85.5626</b> | straight    | 155 mm         |

## Max-Line Arthroscopy Instruments

Shoulder-Instruments, Suture Guide, Suture Hook



Sliding lock handle

Axial Handle

Sutur Retrifer Grasper, serrated



| Sliding lock handle | Axial Handle | Description | Working Length |
|---------------------|--------------|-------------|----------------|
| 50.83.5601          | -            | straight    | 155 mm         |

Sutur Retrifer Grasper, serrated



| Sliding lock handle | Axial Handle | Description                    | Working Length |
|---------------------|--------------|--------------------------------|----------------|
| 50.83.5602          | -            | jaw 45° angled left, straight  | 155 mm         |
| 50.83.5603          | -            | jaw 45° angled right, straight | 155 mm         |

Sutur Retrifer Grasper, serrated



| Sliding lock handle | Axial Handle | Description               | Working Length |
|---------------------|--------------|---------------------------|----------------|
| 50.83.5604          | -            | jaw 15° upbiter, straight | 155 mm         |

King Sutur Retrifer Grasper, cross serrated



| Sliding lock handle | Axial Handle | Description | Working Length |
|---------------------|--------------|-------------|----------------|
| 50.83.5606          | -            | straight    | 155 mm         |

## Max-Line Arthroscopy Instruments

Shoulder-Instruments, Suture Guide, Suture Hook



Handle with positioning

Sliding, Suture Hook



| Handle with positioning | Axial Handle | Description | Working Length |
|-------------------------|--------------|-------------|----------------|
| 50.84.5655              | -            | straight    | 155 mm         |

## Max-Line Arthroscopy Instruments

Hip-Instruments, Punches



Standards handles

Hook Punch, small



| Standard Handle | Description | Tip Profile | Bite Width | Tip Width |
|-----------------|-------------|-------------|------------|-----------|
| 50.80.5701      | straight    | 1.90 mm     | 1.67 mm    | 2.89 mm   |

Hook Punch, small



| Standard Handle | Description       | Tip Profile | Bite Width | Tip Width |
|-----------------|-------------------|-------------|------------|-----------|
| 50.80.5702      | upbiter, straight | 1.90 mm     | 1.67 mm    | 2.89 mm   |

Basket Punch, small



| Standard Handle | Description                       | Tip Profile | Bite Width | Tip Width |
|-----------------|-----------------------------------|-------------|------------|-----------|
| 50.80.5706      | jaw 45° downbiter left, straight  | 1.50 mm     | 1.67 mm    | 2.89 mm   |
| 50.80.5707      | jaw 45° downbiter right, straight | 1.50 mm     | 1.67 mm    | 2.89 mm   |

Basket Punch, small



| Standard Handle | Description | Tip Profile | Bite Width | Tip Width |
|-----------------|-------------|-------------|------------|-----------|
| 50.80.5711      | straight    | 1.93 mm     | 3.20 mm    | 5.00 mm   |
| 50.80.5712      | upswept     | 1.93 mm     | 3.20 mm    | 5.00 mm   |

## Max-Line Arthroscopy Instruments

Hip-Instruments, Punches, Rongeurs



Standards handles

Duckly Punch, medium



| Standard Handle   | Description       | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------------|-------------|------------|-----------|
| <b>50.80.5713</b> | upbiter, straight | 1.93 mm     | 3.20 mm    | 5.00 mm   |



Basket Punch unicorn

| Standard Handle   | Description     | Tip Profile | Bite Width | Tip Width |
|-------------------|-----------------|-------------|------------|-----------|
| <b>50.80.5716</b> | straight, left  | 3.93 mm     | 2.38 mm    | 5.58 mm   |
| <b>50.80.5717</b> | straight, right | 3.93 mm     | 2.38 mm    | 5.58 mm   |

Rongeur Raptor

| Standard Handle   | Description |
|-------------------|-------------|
| <b>50.80.5720</b> | straight    |

## Max-Line Arthroscopy Instruments

Hip-Instruments, Punches, Rongeurs



Rotatable handles



Cigar Handle



### Rotary Basket Punch

| Rotatable Handle  | Cigar Handle      | Description                             | Tip Profile | Bite Width | Tip Width |
|-------------------|-------------------|-----------------------------------------|-------------|------------|-----------|
| <b>50.82.5751</b> | <b>50.81.5751</b> | 2.2 mm, jaw 90° angled left, straight   | 1.60 mm     | 2.20 mm    | 4.10 mm   |
| <b>50.82.5752</b> | <b>50.81.5752</b> | 2.2 mm, jaw 90°, angled right, straight | 1.60 mm     | 2.20 mm    | 4.10 mm   |

## Max-Line Arthroscopy Instruments

Hip-Instruments, Grasping Forceps



Standard handle with ratchet

Grasping Forceps 1x2 tooth

| Standard handle with ratchet | Description |
|------------------------------|-------------|
| 50.83.5760                   | straight    |

Grasping Forceps cupped, serrated



| Standard handle with ratchet | Description      |
|------------------------------|------------------|
| 50.83.5761                   | 2.7 mm, straight |
| 50.83.5762                   | 3.4 mm, straight |



straight

| Standard handle with ratchet | Description |
|------------------------------|-------------|
| 50.83.5763                   | straight    |



## Arthroscopy Instruments - Probes



**50.45.5010**

Hook Probe, pointed



**50.45.5011**

Knife retrograd Ø 5.0 mm, curved 45°



**50.45.5012**

Bajonet Knife



**50.45.5013**

Banana Knife, serrated



**50.45.5014**

Banana Knife, smooth



**50.45.5015**

Spoon Ø 3.0 mm, curved 30°



**50.45.5016**

Ring Curette Ø 3.0 mm, curved 30°



**50.45.5017**

Ring Curette Ø 7.0 mm, curved 30°



**50.45.5018**

Hook Probe Ø 3.0 mm, graduated



**50.45.5019**

Hook Probe Ø 5.0 mm, graduated



**50.45.5020**

Meniscotome Smillie Ø 3.0 mm



**50.45.5021**

Meniscotome Smillie Ø 5.0 mm



## Arthroscopy Instruments - Probes

**50.45.5022**

Meniscotome Smillie Ø 7.0 mm



**50.45.5023**

Scalpel Knife



**50.45.5024**

Rosette Knife Ø 3.0 mm, straight



**50.45.5025**

Rosette Knife Ø 3.0 mm, curved left 30°



**50.45.5026**

Rosette Knife Ø 3.0 mm, curved right 30°



**50.45.5027**

Sickle Knife



**50.45.5028**

Hook Knife Ø 3.0 mm, curved 90°



**50.45.5029**

Rasp narrow straight, fine

**50.45.5030**

Rasp narrow straight, coarse

**50.45.5031**

Oval curette Ø 5.0 mm, curved 15°

**50.45.5032**

Oval curette Ø 3.0 mm, curved 30°

**50.45.5033**

Ring curette Ø 5.0 mm, curved 15°



## Shaver Blades

Multi-use all stainless steel, Non Sterile



Single-use, Non Sterile



Full Radius Resector, Outer and inner blade, non-serrated

| Single-use        | Single-use<br>15° curved | Single-use<br>30° curved | Multi-use all<br>stainless steel | Ø        |
|-------------------|--------------------------|--------------------------|----------------------------------|----------|
| <b>50.25.2501</b> | -                        | -                        | <b>50.35.2501</b>                | Ø 2.5 mm |
| <b>50.25.3501</b> | <b>50.15.3501</b>        | <b>50.30.3501</b>        | <b>50.35.3501</b>                | Ø 3.5 mm |
| <b>50.25.4201</b> | <b>50.15.4201</b>        | <b>50.30.4201</b>        | <b>50.35.4201</b>                | Ø 4.2 mm |
| <b>50.25.5501</b> | <b>50.15.5501</b>        | <b>50.30.5501</b>        | <b>50.35.5501</b>                | Ø 5.5 mm |



Aggressive Full Radius, Outer blade non-serrated, inner blade serrated

| Single-use        | Single-use<br>15° curved | Single-use<br>30° curved | Multi-use all<br>stainless steel | Ø        |
|-------------------|--------------------------|--------------------------|----------------------------------|----------|
| <b>50.25.2502</b> | -                        | -                        | <b>50.35.2502</b>                | Ø 2.5 mm |
| <b>50.25.3502</b> | <b>50.15.3502</b>        | <b>50.30.3502</b>        | <b>50.35.3502</b>                | Ø 3.5 mm |
| <b>50.25.4202</b> | <b>50.15.4202</b>        | <b>50.30.4202</b>        | <b>50.35.4202</b>                | Ø 4.2 mm |
| <b>50.25.5502</b> | <b>50.15.5502</b>        | <b>50.30.5502</b>        | <b>50.35.5502</b>                | Ø 5.5 mm |



Meniscus Cutter, Outer blade serrated, inner blade non-serrated

| Single-use        | Single-use<br>15° curved | Single-use<br>30° curved | Multi-use all<br>stainless steel | Ø        |
|-------------------|--------------------------|--------------------------|----------------------------------|----------|
| <b>50.25.2503</b> | -                        | -                        | <b>50.35.2503</b>                | Ø 2.5 mm |
| <b>50.25.3503</b> | <b>50.15.3503</b>        | <b>50.30.3503</b>        | <b>50.35.3503</b>                | Ø 3.5 mm |
| <b>50.25.4203</b> | <b>50.15.4203</b>        | <b>50.30.4203</b>        | <b>50.35.4203</b>                | Ø 4.2 mm |
| <b>50.25.5503</b> | <b>50.15.5503</b>        | <b>50.30.5503</b>        | <b>50.35.5503</b>                | Ø 5.5 mm |



1 Box: 5 pieces

### Note:

- Maxer Shaver blade single use - sterile packed - minimum purchase 100 pieces.
- Other Shaver blade connection on request.
- Maxer Shaver blades are compatible to Linvatec & Arthrex shaver also.

## Shaver Blades

Trimmer, Outer and inner blade, non-serrated, square window

| Single-use        | Single-use<br>15° curved | Single-use<br>30° curved | Multi-use all<br>stainless steel | Ø        |
|-------------------|--------------------------|--------------------------|----------------------------------|----------|
| <b>50.25.2504</b> | -                        | -                        | <b>50.35.2504</b>                | Ø 2.5 mm |
| <b>50.25.3504</b> | <b>50.15.3504</b>        | <b>50.30.3504</b>        | <b>50.35.3504</b>                | Ø 3.5 mm |
| <b>50.25.4204</b> | <b>50.15.4204</b>        | <b>50.30.4204</b>        | <b>50.35.4204</b>                | Ø 4.2 mm |
| <b>50.25.5504</b> | <b>50.15.5504</b>        | <b>50.30.5504</b>        | <b>50.35.5504</b>                | Ø 5.5 mm |



Aggressive, Meniscus Cutter, Outer and inner blade, serrated

| Single-use        | Single-use<br>15° curved | Single-use<br>30° curved | Multi-use all<br>stainless steel | Ø        |
|-------------------|--------------------------|--------------------------|----------------------------------|----------|
| <b>50.25.2505</b> | -                        | -                        | <b>50.35.2505</b>                | Ø 2.5 mm |
| <b>50.25.3505</b> | <b>50.15.3505</b>        | <b>50.30.3505</b>        | <b>50.35.3505</b>                | Ø 3.5 mm |
| <b>50.25.4205</b> | <b>50.15.4205</b>        | <b>50.30.4205</b>        | <b>50.35.4205</b>                | Ø 4.2 mm |
| <b>50.25.5505</b> | <b>50.15.5505</b>        | <b>50.30.5505</b>        | <b>50.35.5505</b>                | Ø 5.5 mm |



Conical Double Full Radius, Outer and inner blade, non-serrated


| Single-use        | Ø        |
|-------------------|----------|
| <b>50.25.4206</b> | Ø 4.2 mm |



Round Burr

| Single-use        | Multi-use all<br>stainless steel | Ø        |
|-------------------|----------------------------------|----------|
| -                 | <b>50.40.3011</b>                | Ø 3.0 mm |
| <b>50.20.4211</b> | <b>50.40.4211</b>                | Ø 4.2 mm |
| <b>50.20.5511</b> | <b>50.40.5511</b>                | Ø 5.5 mm |





## Shaver Blades

### Tapered Burr

| Single-use        | Multi-use all stainless steel | Ø        |
|-------------------|-------------------------------|----------|
| <b>50.20.5512</b> | <b>50.40.5512</b>             | Ø 5.5 mm |



### Oval Burr

| Single-use        | Multi-use all stainless steel | Ø        |
|-------------------|-------------------------------|----------|
| <b>50.20.5513</b> | <b>50.40.5513</b>             | Ø 5.5 mm |



1 Box : 5 pieces

#### Note:

- Maxer Shaver blade single, use sterile packed, minimum purchase 10 pcs.
- Other Shaver blade connection on request.
- Maxer Shaver blades are compatible to Linvatec & Arthrex shaver also.

## Monopolar Arthroscopy Electrodes

4 mm, Reusable, Autoclavable with Handle



**50.15.1140**

Needle electrode 115 x 1.5 x Ø 0.6 mm, 90° angular, with handle



**50.15.1141**

Needle electrode 115 x 4 x Ø 0.6 mm, 90° angular with handle



**50.15.1142**

Knife electrode 115 x 4 x Ø 1.5 mm, 90° angular with handle



**50.15.1143**

Button electrode with handle, 90° angled, 4 mm, length 115 mm, Ø 1.5 mm



**50.15.1144**

Ball electrode 115 x Ø 3 with handle

## Monopolar Arthroscopy Electrodes

Reusable, Autoclavable for 4 mm pen and 2.36 mm pen, working length 115 mm



|            |            |                                                      |
|------------|------------|------------------------------------------------------|
| 50.15.2151 | 50.15.3151 | Needle electrode 90° angular,<br>Ø 1.5 mm x Ø 0.6 mm |
|------------|------------|------------------------------------------------------|



|            |            |                                                      |
|------------|------------|------------------------------------------------------|
| 50.15.2152 | 50.15.3152 | Needle electrode 90° angular,<br>Ø 4.0 mm x Ø 0.6 mm |
|------------|------------|------------------------------------------------------|



|            |            |                                                      |
|------------|------------|------------------------------------------------------|
| 50.15.2153 | 50.15.3153 | Needle electrode 90° angular,<br>Ø 3.0 mm x Ø 1.5 mm |
|------------|------------|------------------------------------------------------|



|            |            |                             |
|------------|------------|-----------------------------|
| 50.15.2154 | 50.15.3154 | Brush electrode angular 70° |
|------------|------------|-----------------------------|



|            |            |                                      |
|------------|------------|--------------------------------------|
| 50.15.2155 | 50.15.3155 | Button electrode, Ø 3.0 mm, straight |
|------------|------------|--------------------------------------|



|            |            |                                    |
|------------|------------|------------------------------------|
| 50.15.2156 | 50.15.3156 | Ball electrode hollow, angular 70° |
|------------|------------|------------------------------------|



|            |            |                                  |
|------------|------------|----------------------------------|
| 50.15.2157 | 50.15.3157 | Ball electrode flat, angular 70° |
|------------|------------|----------------------------------|



## Monopolar - High Frequency Cable



|                   |                                                          |
|-------------------|----------------------------------------------------------|
| <b>20.55.0003</b> | High frequency cable, monopolar, male / male, 3 m length |
|-------------------|----------------------------------------------------------|

|                   |                                                          |
|-------------------|----------------------------------------------------------|
| <b>20.55.0005</b> | High frequency cable, monopolar, male / male, 5 m length |
|-------------------|----------------------------------------------------------|

|                   |                                                            |
|-------------------|------------------------------------------------------------|
| <b>20.55.1003</b> | High frequency cable, monopolar, male / female, 3 m length |
|-------------------|------------------------------------------------------------|

|                   |                                                            |
|-------------------|------------------------------------------------------------|
| <b>20.55.1005</b> | High frequency cable, monopolar, male / female, 5 m length |
|-------------------|------------------------------------------------------------|


|                   |                                                              |
|-------------------|--------------------------------------------------------------|
| <b>20.55.2003</b> | High frequency cable, monopolar, female / female, 3 m length |
|-------------------|--------------------------------------------------------------|

|                   |                                                                                     |
|-------------------|-------------------------------------------------------------------------------------|
| <b>20.55.3003</b> | High frequency cable, monopolar, female / male = "Big Jack" (Valleylab), 3 m length |
|-------------------|-------------------------------------------------------------------------------------|



## Bipolar Arthroscopy Electrodes

Reusable, Autoclavable for use in SALINE SOLUTION



|            |            |                                                         |
|------------|------------|---------------------------------------------------------|
| 50.15.1701 | 50.15.1151 | Bipolar vaporization electrode, "Phazer"<br>angular 90° |
|------------|------------|---------------------------------------------------------|



|            |            |                                                         |
|------------|------------|---------------------------------------------------------|
| 50.15.1702 | 50.15.1152 | Bipolar vaporization electrode, "Phazer"<br>angular 70° |
|------------|------------|---------------------------------------------------------|



|            |            |                                                               |
|------------|------------|---------------------------------------------------------------|
| 50.15.1703 | 50.15.1153 | Bipolar needle electrode,<br>Ø 2.0 mm X Ø 0.6 mm, angular 90° |
|------------|------------|---------------------------------------------------------------|



|            |            |                                                            |
|------------|------------|------------------------------------------------------------|
| 50.15.1704 | 50.15.1154 | Bipolar needle electrode,<br>Ø 4.0 X Ø 0.6 mm, angular 90° |
|------------|------------|------------------------------------------------------------|



|            |            |                           |
|------------|------------|---------------------------|
| 50.15.1705 | 50.15.1155 | Bipolar electrode "Cobra" |
|------------|------------|---------------------------|



## Bipolar Instrument Cable

Cable for instruments, 3 m length

|            | Bipolar-Units                                                                                                   | Instruments                                                                                        |
|------------|-----------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| 20.67.2003 | <br>Erbe/Wisap/Storz           | <br>flat plug   |
| 20.67.2103 | <br>GIMMI / Martin / Berchtold | <br>flat plug   |
| 20.67.2203 | <br>Valleylab                | <br>flat plug |



## NOTES

[illegible]



UNITS

LAPAROSCOPY

UROLOGY

HYSTEROSCOPY

ARTHROSCOPY

ENT

PEDIATRIC

SPINE ENDOSCOPY

Distributed by:

**Maxer Medizintechnik GmbH**

Untere Hauptstr. 34/1

78573 Wurmlingen

Germany

Registered Office: Vogesenstr. 17

78549 Spaichingen

Germany

Email: [info@maxerendoscopy.com](mailto:info@maxerendoscopy.com)

Website: [www.maxerendoscopy.com](http://www.maxerendoscopy.com)