

Maxer
GERMANY

LAPAROSCOPY



M LAPAROSCOPY
SYSTEM

maxerendoscopy.com



ISO Certificate

Quality management system as per
ISO 13485 :2012 Medical Device Certification
for development & Production & Endoscopic
instrument & Equipments

CE Certificate

As per council directive 93/42/EEC
Concerning medical devices



About Maxer

Maxer Medizintechnik GmbH has established itself as a symbol of creativity and latest technology and is an internationally qualified supplier of world class medical instruments and devices.

Specialized in the field of Minimal Invasive Surgery, Maxer manufactures wide range of equipments such as Laparoscopy, Hysteroscopy, Urology, Pediatric Surgery, Arthroscopy, Spine Endoscopy, ENT, etc.

Maxer has been one of the foremost companies to bring out HD technology in MIS field. In 2004, Maxer brought out first endoscopy camera with Matrix Metering. In 2008, Maxer brought out HD Telescopes. In 2009, Maxer launched first Networked Endoscopy Camera System. In 2010, Maxer introduced pathbreaking 16 Fr Resectoscope for Office Hysteroscopy. In 2011, Maxer launched LED Headlight with 200,000 Lux illumination, again a state of art in this segment. In 2012, Maxer launched LED Light sources to keep Maxer customers always at the cutting edge of technology. Maxer has also launched HD cameras with integrated Full HD recording. 3 CMOS HD cameras with 6 megapixel camera head will be launched soon.

Owners of Maxer have been working in Medical market for more than 20 years and have installed thousands of sets with high level of customer satisfaction. This background has given Maxer profound credibility in the field of endoscopy.

Due to vibrant co-operation with end users, Maxer is able to continuously improve upon the product range and also able to meet special modification demands. Maxer believes that growth and innovation go hand-in-hand.









Maxer sells complete system with a focus on product quality and customers satisfaction. Maxer provides sales and service training to it's distributor in order to be successful.

Maxer believes that sound ethics with professionalism will only give success to the company. Maxer's growth and success are based on long term partnerships with customers, distributors. Maxer currently exports it's products to Europe, Asia, Africa.

Maxer is headed by Mr. Giovanni Miranda in the capacity as a Managing Director.

How to reach to Maxer - Maxer is located in Würlingen (sales office) within the medical instruments manufacturing region around Tuttlingen in South Germany. Maxer also has registered office at Spaichingen.




INDEX

Laparoscopes & Information	L - 4 to 7
Trocars & Accessories	L - 8 to 20
Veress (Insufflation) Needles	L - 21
Suction & Irrigation Cannulas	L - 22 to 24
Cholangiography Fixation Forceps	L - 25
Claw Forceps	L - 26
Clip Applying Forceps	L - 27 to 29
Endo Retractor & Stand holder for Liver Retractor	L - 30 to 31
Needle Holder	L - 32 to 33
Knot Guide	L - 34
Ring Applicator	L - 35
Uterus Manipulator	L - 36 to 38
Myoma Drill & Uterin Probe	L - 39
Suction & Irrigation Handles, Electrodes	L - 40
Monopolar Electrodes	L - 41 to 43
Monopolar Cable	L - 44
Cleaning Instruments	L - 45
Bipolar Electrodes	L - 46 to 50
Bipolar Cable	L - 51
SLIMLINE Bipolar	L - 52
SLIMLINE Bipolar Cables	L - 53
Bipolar Coagulation Electrodes Cable	L - 54 to 55
Bariatric Laparoscopy	L - 56
Graspers, Dissectors & Scissors	L - 57 to 70

Laparoscopes from Maxer, Germany


Maxer has been setting innovative standards for many years. Maxer produces laser welded endoscopes which can withstand numerous autoclaving cycles. Its glass cone technology is also unique. Maxer endoscopes are known worldwide for third generation optics which give high resolution and brightness all over the image. New High Definition Laparoscopes are developed with better optical clarity and depth of focus.



Ø 5.0 mm	0°	30°	45°
Ø 5.5 mm			
Ø 10.0 mm	0°	30°	45°


Sapphire tip -

For 3 times more scratch resistance,
safe with shaver blades



Laser welded joints -

No chance of leakage



Glass cone -

For 5 mm & smaller scopes
For longer life of fibers

Complete stainless steel -
Increased sturdiness



- High Definition Laparoscopes have better optics to give enhanced clarity & brightness
- High-quality, multi coated (anti-reflective) optical system
- Highly refractive glass for optimal color rendering even when the image is magnified
- Extremely precise centering of the optical and mechanical components for an image that is brilliant and high contrast even around the edges
- Triple tube design for higher stability with many types of endoscopes
- Nitrogen filled optical system to avoid interior condensation
- High resolution scopes, Sapphire cover glass on both ends,
- Glass cone to collect light from F. O. Cable, Autoclavable

Laparoscopes

Ø 5 mm & Ø 5.5 mm



HM version 5 mm

0°	30°	45°	Ø mm	Working Length
15.10.5100	15.10.5130	15.10.5145	5.5	300 mm

HM version 5.5 mm

0°	30°	45°	Ø mm	Working Length
15.10.5600	15.10.5630	15.10.5645	5.5	300 mm

HM HD version 5.5 mm

0°	30°	45°	Ø mm	Working Length
15.10.6600	15.10.6630	15.10.6645	5.5	300 mm

HD bariatric laparoscope version 5.5 mm

0°	30°	45°	Ø mm	Working Length
15.10.2500	15.10.2530	15.10.2545	5.5	425 mm

HM = High Magnification offers magnified image for the same working area

HD = High Definition scopes offer better optical clarity, brightness & depth of focus

Laparoscopes

Ø 10 mm



Standard Version 10 mm

0°	30°	45°	Ø mm	Working Length
15.10.1000	15.10.1030	15.10.1045	10.0	330 mm

HD Version 10 mm

0°	30°	45°	Ø mm	Working Length
15.10.1100	15.10.1130	15.10.1145	10.0	330 mm

HD⁺ Laparoscopes



These are newly developed telescopes which provide outstanding sharpness at periphery as well as at center of the image. These optics are designed to match latest range of Full HD Cameras from Maxer.

HD⁺ Version 10 mm

0°	30°	45°	Ø mm	Working Length
15.10.1101	15.10.1131	15.10.1146	10.0	330 mm

Bariatric Laparoscope Version 10 mm

0°	30°	45°	Ø mm	Working Length
15.10.2100	15.10.2130	15.10.2145	10.0	450 mm

Laparoscopes

Ø 10 mm & Ø 11 mm



Standard Version


90° angled	Ø mm	Working Channel	Working Length
15.10.9000	10.0	Ø 5.1 mm	270 mm
15.10.9011	11.0	Ø 6.1 mm	270 mm



Standard Version

45° angled	Ø mm	Working Channel	Working Length
15.10.4500	10.0	Ø 5.1 mm	270 mm
15.10.4511	11.0	Ø 6.1 mm	270 mm


Trocars and accessories



Automatic trocar sleeve with stop-cock for insufflation & trocar mandrin with pyramidal tip

Sleeve with trocar, compl	Trocar sleeve only	Trocar only	Ø mm	Working Length
20.13.1035	20.13.0035	20.12.1035	3.5	70 mm
20.13.1050	20.13.0050	20.12.1050	5.0	92 mm
20.13.1055	20.13.0055	20.12.1055	5.5	95 mm
20.13.1060	20.13.0060	20.12.1060	6.0	95 mm
20.13.1070	20.13.0070	20.12.1070	7.0	95 mm
20.13.1080	20.13.0080	20.12.1080	8.0	105 mm
20.13.1100	20.13.0100	20.12.1100	10.0	97 mm
20.13.1110	20.13.0110	20.12.1110	11.0	101 mm
20.13.1125	20.13.0125	20.12.1125	12.5	98 mm


Trocars and accessories



Automatic Trocar sleeve with stop cock for insufflation and trocar mandrin with conical tip

Sleeve with trocar compl	Trocar sleeve only	Trocar only	Ø mm	Working Length
20.13.2035	20.13.0035	20.12.2035	3.5	70 mm
20.13.2050	20.13.0050	20.12.2050	5.0	92 mm
20.13.2055	20.13.0055	20.12.2055	5.5	95 mm
20.13.2060	20.13.0060	20.12.2060	6.0	95 mm
20.13.2070	20.13.0070	20.12.2070	7.0	95 mm
20.13.2080	20.13.0080	20.12.2080	8.0	105 mm
20.13.2100	20.13.0100	20.12.2100	10.0	97 mm
20.13.2110	20.13.0110	20.12.2110	11.0	101 mm
20.13.2125	20.13.0125	20.12.2125	12.5	98 mm


Trocars and accessories



Multiple function trocar sleeve with manual opening lever, stop cock for insufflation & trocar mandrin with pyramidal tip

Sleeve with trocar compl	Trocar sleeve only	Trocar only	Ø mm	Working Length
20.14.1035	20.14.0035	20.12.1035	3.5	70 mm
20.14.1050	20.14.0050	20.12.1050	5.0	92 mm
20.14.1055	20.14.0055	20.12.1055	5.5	95 mm
20.14.1060	20.14.0060	20.12.1060	6.0	95 mm
20.14.1070	20.14.0070	20.12.1070	7.0	95 mm
20.14.1080	20.14.0080	20.12.1080	8.0	105 mm
20.14.1100	20.14.0100	20.12.1100	10.0	97 mm
20.14.1110	20.14.0110	20.12.1110	11.0	101 mm
20.14.1125	20.14.0125	20.12.1125	12.5	98 mm
20.14.1130	20.14.0130	20.12.1130	13.0	97 mm
20.14.1150	20.14.0150	20.12.1150	15.0	95 mm
20.14.1200	20.14.0200	20.12.1200	20.0	95 mm

Trocars and accessories




Multiple function Trocar sleeve with manual opening lever, stop cock for insufflation and trocar mandrin with conical tip

Sleeve with trocar compl	Trocar sleeve only	Trocar only	Ø mm	Working Length
20.14.2035	20.14.0035	20.12.2035	3.5	70 mm
20.14.2050	20.14.0050	20.12.2050	5.0	92 mm
20.14.2055	20.14.0055	20.12.2055	5.5	95 mm
20.14.2060	20.14.0060	20.12.2060	6.0	95 mm
20.14.2070	20.14.0070	20.12.2070	7.0	95 mm
20.14.2080	20.14.0080	20.12.2080	8.0	105 mm
20.14.2100	20.14.0100	20.12.2100	10.0	97 mm
20.14.2110	20.14.0110	20.12.2110	11.0	101 mm
20.14.2125	20.14.0125	20.12.2125	12.5	98 mm
20.14.2130	20.14.0130	20.12.2130	13.0	97 mm
20.14.2150	20.14.0150	20.12.2150	15.0	95 mm
20.14.2200	20.14.0200	20.12.2200	20.0	95 mm

Trocars and accessories

Hasson Trocars



Hasson automatic trocar with suture holder fixed at cone, stop-cock for insufflation & trocar mandrin with blunt tip

HASSON trocar compl.	Sleeve with cone, only	Trocar blunt, only	Ø mm
20.15.3055	20.15.0055	20.12.3055	5.5
20.15.3110	20.15.0110	20.12.3110	11.0
20.15.3112	20.15.0125	20.12.3125	12.5

Fixable Cone only, for HASSON trocar, with suture holder, fit for all trocar sleeves

Art No.	Ø mm
20.15.4035	3.5
20.15.4055	5.5
20.15.4010	10.0
20.15.4011	11.0
20.15.4125	12.5


Fixable Cone with thread only, for HASSON trocar, fit for all trocar sleeves, without suture holder

20.15.5110	11.0
-------------------	------



Trocars and accessories

Thoracic Trocar



Rigid Thoracic Trocars

Sleeve with trocar compl	Trocar sleeve only	Trocar only	Ø mm	Working Length
20.16.5060	20.16.0060	20.12.5060	6.0	50 mm
20.16.5110	20.16.0110	20.12.5110	11.0	50 mm
20.16.5130	20.16.0130	20.12.5130	13.0	50 mm



Flexible thoracic trocars

Sleeve with trocar compl	Trocar sleeve only	Trocar only	Ø mm	Working Length
20.16.6060	20.16.7060	20.12.6060	6.0	100 mm
20.16.6061	20.16.7061	20.12.6061	6.0	70 mm
20.16.6110	20.16.7110	20.12.6110	11.0	100 mm
20.16.6111	20.16.7111	20.12.6111	11.0	70 mm

Flexible plastic tubes box of 5 each

Art No.	for trocar sleeve Ø
20.16.8061	6.0 / 100
20.16.8062	6.0 / 70
20.16.8110	11.0 / 100
20.16.8111	11.0 / 70

Trocars and accessories



Security trocar, pyramidal tip (detachable) useable with all trocar sleeves

Art No.	Ø mm	Working Length
20.12.8050	5.0	157 mm
20.12.8055	5.5	160 mm
20.12.8080	8.0	171 mm
20.12.8100	10.0	167 mm
20.12.8110	11.0	171 mm
20.12.8125	12.5	171 mm



Security trocar, oblique (detachable) useable with all trocar sleeves

20.12.9050	5.0	157 mm
20.12.9055	5.5	160 mm
20.12.9080	8.0	171 mm
20.12.9100	10.0	167 mm
20.12.9110	11.0	171 mm
20.12.9125	12.5	171 mm

Trocars and accessories



Security trocar, oblique useable with all trocar sleeves

Art No.	Ø mm	Working Length
20.17.9050	5.0	157 mm
20.17.9055	5.5	160 mm
20.17.9080	8.0	171 mm
20.17.9100	10.0	167 mm
20.17.9110	11.0	171 mm
20.17.9125	12.5	171 mm



Security Trocar, automatic locking, obturator tip "Knife Design"

Article No.	Ø mm	Working Length
20.10.4055	5.5	160 mm
20.10.4110	11	171 mm
20.10.4125	12.5	171 mm


Trocars and accessories

Trocar tubes, interchangeable, for automatic trocars

		For Ø mm trocar sleeve	Working Length
20.13.4035	-	3.5 mm	70 mm
20.13.4050	20.13.5050	5.0 mm	92 mm
20.13.4055	20.13.5055	5.5 mm	95 mm
20.13.4060	20.13.5060	6.0 mm	95 mm
20.13.4080	20.13.5080	8.0 mm	105 mm
20.13.4100	20.13.5100	10.0 mm	97 mm
20.13.4110	20.13.5110	11.0 mm	101 mm
20.13.4125	20.13.5125	12.5 mm	98 mm
20.13.4130	-	13.0 mm	97 mm
20.13.4150	-	15.0 mm	95 mm
20.13.4200	-	20.0 mm	95 mm

Trocar stop for trocar sleeve

Length of Thread 35 mm	For Ø mm Trocar Sleeve
20.14.3505	5.0 mm
20.14.3555	5.5 mm
20.14.3510	10.0 mm
20.14.3511	11.0 mm
20.14.3512	12.5 mm



Trocars and accessories



Reducer sleeve

Art No.	Trocar sleeve Ø mm	Instrument Ø mm	Working Length
20.18.0504	5.0	4.0	175 mm
20.18.1055	10.0	5.5	190 mm
20.18.1155	11.0	5.5	190 mm
20.18.1255	12.5	5.5	190 mm
20.18.1210	12.5	10.0	190 mm
20.18.1555	15.0	5.5	190 mm
20.18.1510	15.0	10.0	190 mm
20.18.2055	20.0	5.5	190 mm
20.18.2010	20.0	10.0	190 mm
20.18.2015	20.0	15.0	190 mm

Trocars and accessories


Reducer, removable, with Holder

Art No.	from Ø mm	to Ø mm
20.19.0855	8.0	5.5
20.19.1055	10.0	5.5
20.19.1155	11.0	5.5
20.19.1255	12.5	5.5
20.19.1210	12.5	10.0



Reducer, removable

Art No.	from Ø mm	to Ø mm
20.20.1055	10.0	5.5
20.20.1155	11.0	5.5
20.20.1255	12.5	5.5
20.20.1210	12.5	10.5
20.20.1555	15.0	5.5
20.20.1510	15.0	10.0
20.20.2055	20.0	5.5
20.20.2010	20.0	10.0



Trocars and accessories



Dilation - Set, consisting of dilator sleeve and dilator

20.42.0510	Ø 5.0 mm	Ø 10.0 mm
20.42.0511	Ø 5.0 mm	Ø 11.0 mm
20.42.1015	Ø 10.0 mm	Ø 15.0 mm
20.42.0515	Ø 5.0 mm	Ø 15.0 mm
20.42.0520	Ø 5.0 mm	Ø 20.0 mm
20.42.1020	Ø 10.0 mm	Ø 20.0 mm



Gall Stone collector, Ø 20 mm

20.42.2055	For instruments Ø 5.5 mm
20.42.2010	For instruments Ø 10 mm



20.42.3000	Retractor for gall bladder extraction with 3 blades, jaw length 60 mm
-------------------	---

Trocars and accessories



Outer



Inner

Sealing Caps / Outer & Inner Washers for Trocars

20.90.1001	Trocar washer outer Ø 11 mm
20.90.1002	Trocar washer inner Ø 11 mm
20.90.1003	Trocar washer outer Ø 5.5 mm
20.90.1004	Trocar washer inner Ø 5.5 mm
20.90.1013	Trocar washer outer Ø 6 mm
20.90.1014	Trocar washer inner Ø 6 mm
20.90.1005	Trocar washer outer Ø 7 mm
20.90.1006	Trocar washer inner Ø 7 mm
20.90.1007	Trocar washer outer Ø 12 mm
20.90.1008	Trocar washer Ø 12 mm
20.90.1009	Trocar washer outer Ø 15 mm
20.90.1010	Trocar washer inner Ø 15 mm

O-Rings

20.90.4001	For trocar Ø 5.5 & 6 mm
20.90.4002	For trocar Ø 7 & 8 mm
20.90.4003	For trocar Ø 10 & 11 mm
20.90.4004	For trocar Ø 12, 12.5 & 13 mm
20.90.4005	For trocar Ø 15 mm

Veress (Insufflation) Needles



20.10.1055	Veress needle Ø 2.0 mm, WL 55 mm
20.10.1070	Veress needle Ø 2.0 mm, WL 70 mm
20.10.1080	Veress needle Ø 2.0 mm, WL 80 mm
20.10.1100	Veress needle Ø 2.0 mm, WL 100 mm
20.10.1120	Veress needle Ø 2.0 mm, WL 120 mm
20.10.1130	Veress needle Ø 2.0 mm, WL 130 mm
20.10.1140	Veress needle Ø 2.0 mm, WL 140 mm
20.10.1150	Veress needle Ø 2.0 mm, WL 150 mm
20.10.1180	Veress needle Ø 2.0 mm, WL 180 mm
20.10.1200	Veress needle Ø 2.0 mm, WL 200 mm



20.10.2100	Veress needle high flow, Ø 2.7 mm, WL 100 mm
20.10.2120	Veress needle high flow, Ø 2.7 mm, WL 120 mm
20.10.2150	Veress needle high flow, Ø 2.7 mm, WL 150 mm
20.10.2170	Veress needle high flow, Ø 2.7 mm, WL 170 mm

Suction & Irrigation Cannulas



20.21.2001

Suction & Irrigation Handle Ø 5 mm, with sliding valves



20.21.2005

Tube, interchangeable, Ø 5 mm



20.21.2010

Tube, interchangeable, Ø 10 mm

20.21.2000

Suction & Irrigation set (handle 20.21.2001 & interchangeable tubes Ø 5 mm 20.21.2005 & Ø 10 mm 20.21.2010)



20.21.3102

Suction & Irrigation tube Ø 5 mm with two trumpet valves graduated, channel for laser and square handle, WL 34 cm




20.21.3103

Laser guide with Luer-Lock connection, WL 39 cm, for handle 20.21.3102

Suction & Irrigation Cannulas


20.21.3001

Handle with host join and two trumpet valves



20.21.3002

Handle with LL and two trumpet valves



20.21.3003

Handle with sliding valves



20.21.3330

Suction tube Ø 5 mm with 4 holes, Standard, L 330 mm



20.21.3331

Suction tube Ø 5 mm with 16 holes, blunt, L 327.5 mm



20.21.3332

Suction tube Ø 10 mm with 32 holes, blunt, L 330 mm



Suction & Irrigation Cannulas



20.21.3333

Suction tube Ø 5 mm with 48 holes, L 330 mm



20.21.3334

Suction tube Ø 5 mm with 48 holes, blunt, L 500 mm



20.21.3335

Aspiration needle Ø 5 mm, L 330 mm



20.21.3336

Aqua dissector cannula Ø 5 mm, L 330 mm

SET of Suction & Irrigation tubes, consisting of : 7 changeable tubes and 1 handle

20.21.4001

Set with handle **20.21.3001**

20.21.4002

Set with handle **20.21.3002**

20.21.4003

Set with handle **20.21.3003**

Cholangiography Fixation Forceps



20.26.0533

Cholangiography fixation forceps Ø 5 mm, with inside ratchet, WL 33 cm

20.26.0528

Cholangiography fixation forceps Ø 5 mm, with inside ratchet, WL 28 cm



20.26.1533

Cholangiography Fixation Claw Forceps Ø 5 mm, handle with retaining spring, WL 33 cm

Claw Forceps

Scissors, Spoons with Rinsing Channel



20.30.1033

Claw grasping forceps Ø 10 mm, handle with retaining spring, WL 33 cm



20.30.2033

Scissors one blade toothed, Ø 10, handle without retaining spring, WL 33 cm



20.30.3033

Spoon forceps Ø 10 mm, handle with retaining spring, WL 33 cm


Clip Applying Forceps



Clip applying forceps Ø 10 mm, WL 33 cm, 360° rotatable, with LL rinsing channel, for TITANIUM Clips

Art No.	For Clips - Type
20.38.1031	Ethicon - medium
20.38.1032	Ehticon - medium - large
20.38.1033	Ethicon - large
20.38.1034	Pilling/Weck - medium
20.38.1035	Pilling/Weck - medium-large
20.38.1036	SLS - medium
20.38.1037	SLS - medium - large
20.38.1038	SLS - large
20.38.1039	HORIZON - medium
20.38.1040	HORIZON - medium-large
20.38.1041	HORIZON - large

Clip Applying Forceps



Clip applying forcep Ø 10 mm, WL 33 cm, detachable, 360° rotatable, for TITANIUM Clips

Forceps complete Art No.	Insert only Art No.	For Clips - Type
20.38.2031	20.38.3031	Ethicon - medium
20.38.2032	20.38.3032	Ehticon - medium-large
20.38.2033	20.38.3033	Ethicon - large
20.38.2034	20.38.3034	Pilling/Weck - medium
20.38.2035	20.38.3035	Pilling/ Weck - medium - large
20.38.2036	20.38.3036	SLS - medium
20.38.2037	20.38.3037	SLS - medium - large
20.38.2038	20.38.3038	SLS - large
20.38.2039	20.38.3039	HORIZON - medium
20.38.2040	20.38.3040	HORIZON - medium
20.38.2041	20.38.3041	HORIZON - large

Clip Applying Forceps

Ø 5 mm Rotatable



Clip applying forceps Ø 5 mm, WL 25 cm and 33 cm, detachable and 360° rotatable for TITANIUM Clips

Forceps, complete 25 cm	Forceps, complete 33 cm	For Clips - Type
20.38.4031	20.38.5031	Ethicon - small
20.38.4032	20.38.5032	Ethicon - medium
20.38.4033	20.38.5033	Pilling/Weck - small
20.38.4034	20.38.5034	Pilling/Weck - medium
20.38.4035	20.38.5035	SLS - small
20.38.4036	20.38.5036	SLS - small - medium
20.38.4037	20.38.5037	SLS - medium
20.38.4038	20.38.5038	HORIZON - small
20.38.4039	20.38.5039	HORIZON - medium

Endo Retractors

with irrigation device



20.40.5134

Endo Retractor 3 fingers Ø 5 mm, rigid, WL 34 cm



20.40.5138

Endo Retractor 5 fingers Ø 10 mm, rigid, WL 38 cm



20.40.5234


Endo Retractor 3 fingers Ø 5 mm, angleable, WL 34 cm



20.40.5238

Endo Retractor 5 fingers Ø 10 mm, angleable, WL 38 cm

Endo Retractors



20.40.5334

Flexible Liver Retractor for Ø 5 mm trocar

Can be folded using outside knob.

The retractor offers advantage of superior strength and the rounded edges cause minimum damage to the Liver



12.03.0010

Stand holder for Liver Retractor

12.03.0011

Table clamp

Needle Holders

2 branch handle & LL rinsing channel



20.35.6033

Endo Needle holder Ø 5.0 mm, straight, WL 33 cm



20.35.6133

Endo Needle holder Ø 5.0 mm, straight plain, WL 33 cm



20.35.6233

Endo Needle holder Ø 5.0 mm, straight nose, WL 33 cm



20.35.6333

Endo Needle holder Ø 5.0 mm, right curved, WL 33 cm



20.35.6433

Endo Needle holder Ø 5.0 mm, left curved, WL 33 cm



Needle Holders

2 branch handle, LL rinsing channel



20.35.6666

Endo Needle Holder tungsten carbide Ø 5 mm, straight, with overload protection, L 33 cm



20.35.6833

Endo Needle Holder tungsten carbide Ø 5 mm, conical straight, light weight, L 33 cm

20.35.6933

Endo Needle Holder tungsten carbide Ø 5 mm, right curved, light weight, L 33 cm

20.35.7033

Endo Needle Holder tungsten carbide Ø 5 mm, left curved, light weight, L 33 cm



20.28.3400

Debakey grasping forcep double action, tip 2.8 mm, Ø 5 mm, pencil handle, metallic tubing, luer lock connection, WL 34 cm

Knot Guide



20.37.5133

Knot guide with slot Ø 4.5 mm, WL 33 cm



20.37.5203

Knot guide with eye Ø 4 mm, open, WL 33 cm

20.37.5205

Knot guide with eye Ø 4 mm, open, WL 45 cm



20.37.6203

Knot guide with eye Ø 4 mm, closed, WL 33 cm

20.37.6205

Knot guide with eye Ø 4 mm, closed, WL 45 cm



20.37.7203

Knot guide with eye Ø 3 mm, curved open, WL 33 cm

20.37.7205

Knot guide with eye Ø 5 mm, curved open, WL 33 cm



Ring Applicators



20.75.0601

Ring applicator for tube sterilisation, outer Ø 6 mm for single puncture for use with operating laparoscopes with 6 mm instrument channels



20.75.0602


Ring applicator for tube sterilisation, outer Ø 6 mm for second puncture for use with trocar sleeves 7.0 mm



20.35.6043

Port Closure Forceps

Uterus Manipulator



20.70.0040

Uterus Manipulator only

20.70.1040

Adaptor Ø 8 mm x 80 mm



20.70.2040

Adaptor Ø 10 mm x 100 mm



20.70.3040

Adaptor Ø 3 mm x 45 mm



20.70.4040

Adaptor Ø 3 mm x 55 mm



20.70.5040

Adaptor Cone Short



20.70.6040


Adaptor Cone Long



20.70.0040S

Uterus Manipulator complete set (Manipulator + 6 adaptors)

Uterus Manipulator



The Uterus Manipulator was developed for the manipulation and presentation of the uterus during a complete laparoscopic hysterectomy, as well as protecting the bladder, ureter, and rectum during the electrosurgical resection of the uterus. The reusable instrument is fully modular and can fit all anatomical conditions as vaginal delineating tubes and intrauterine elements of various sizes are available. Once the instrument is in place the uterus can be moved in all directions, thus providing optimal exposure of all anterior, posterior and lateral aspects of the uterus.

The vaginal delineating tube is slipped over the instrument shaft in a forward direction and thus allows the separation of the vagina whilst offering protection to the bladder, ureter and rectum. The electrosurgical resection can be easily and safely accomplished around the ceramic ring. The ceramic ring does not deteriorate when electrosurgical energy is applied and thus can be reused many times.

The vaginal delineating tube is available in 3 different sizes to fit the specific anatomical conditions. The assembled instrument is comparatively light in weight and well balanced – this reduces fatigue during use.

With the tenaculum forceps tightly fixed to the cervix, the uterus can now easily be moved for presentation and preparation purposes.

The complete instrument set consists of several components which can be easily dismantled, cleaned, sterilised and reassembled. The high-grade materials are very durable and sterilisation is trouble-free.

Uterus Manipulator



20.70.0041	Universal handle, rotatable
20.70.0141	Shaft with fixation nut
20.70.0142	Inner shaft with locking pin
20.70.0143	Sliding and fixation shaft
20.70.0144	Cervical tenaculum

20.70.1041	Intrauterine element 4 x 40 mm
20.70.1042	Intrauterine element 5 x 50 mm
20.70.1043	Intrauterine element 6 x 60 mm
20.70.1044	Intrauterine element 7 x 70 mm
20.70.1045	Intrauterine element 8 x 80 mm

20.70.3500	Tube 35 mm with insulating ring
20.70.4000	Tube 40 mm with insulating ring
20.70.4500	Tube 45 mm with insulating ring



A complete set consists of 1 of each of the above items.

Myoma Drill & Uterin Probe



20.71.0533

Myoma Drill, Ø 5 mm, WL 33 cm

20.71.1033

Myoma Drill, Ø 10 mm, WL 33 cm



20.72.0020

COHEN Intra Uterin Probe with 1 small cone and 1 large cone



20.73.0050

Tube dialator Ø 5 mm, WL 34 cm

Suction & Irrigation Handles, Electrodes



20.45.2002

Pistol handle with 2 way valve for adaptable HF Electrodes with Suction & Irrigation channel, rotating adaptor with lateral HF Connector and 2 universal hose connectors

Adaptable monopolar HF-Electrodes

with Suction & Irrigation channel, Ø 5 mm, WL 340 mm, for Handle **20.45.2002**



20.50.1001

Hook electrode thick



20.50.1002

Right angle electrode



20.50.1003

Spatula electrode



20.50.1004

Ballpoint electrode



20.50.1007

Needle electrode



Monopolar Electrodes



HF-Electrodes, monopolar, with HF connector and Straight suction / flush port, Ø 5 mm

WL 33 cm	WL 45 cm	
----------	----------	--

20.50.3001	20.50.3101	Hook electrode
-------------------	-------------------	----------------



20.50.3002	20.50.3102	Right angle electrode
-------------------	-------------------	-----------------------



20.50.3003	20.50.3103	Spatula electrode
-------------------	-------------------	-------------------



20.50.3004	20.50.3104	Ballpoint electrode
-------------------	-------------------	---------------------



20.50.3007	20.50.3107	Needle electrode
-------------------	-------------------	------------------



Monopolar Electrodes



20.50.4101

Spatula electrode, Ø 5 mm, WL 33 cm



20.50.4102

Ballpoint electrode, Ø 5 mm, WL 33 cm



20.50.4103

Hook electrode "J", Ø 5 mm, WL 33 cm



20.50.4104

Hook electrode "L", Ø 5 mm, WL 33 cm



20.50.4105

Needle electrode, Ø 5 mm, WL 33 cm



20.50.4100

Comfort handle for electrodes



Monopolar Electrodes

With trumpet valve handle (Not detachable)



20.50.5001	Spatula electrode monopolar Ø 5 mm, WL 33 cm
-------------------	--



20.50.5002	suction irrigation tube Ø 5 mm, WL 33 cm
-------------------	--



20.50.5003	J Hook electrode monopolar Ø 5 mm, WL 33 cm
-------------------	---



20.50.5004	L hook electrode monopolar Ø 5 mm, WL 33 cm
-------------------	---



20.50.5005	Aspiration needle suction tube Ø 5 mm, WL 33 cm
-------------------	---



Monopolar Ceramic Electrodes



20.50.2001	Spatula electrode, monopolar Ø 5 mm, WL 33 cm
-------------------	---




20.50.2003	Hook electrode "J", monopolar Ø 5 mm, WL 33 cm
-------------------	--



20.50.2004	Hook electrode "L", monopolar Ø 5 mm, WL 33 cm
-------------------	--



High Frequency Cable, monopolar



20.55.0003	High frequency cable, monopolar, male / male, 3 m length
20.55.0005	High frequency cable, monopolar, male / male, 5 m length
20.55.1003	High frequency cable, monopolar, male / female, 3 m length
20.55.1005	High frequency cable, monopolar, male / female, 5 m length
20.55.2003	High frequency cable, monopolar, female / female, 3 m length
20.55.3003	High frequency cable, monopolar, female / male = "Big Jack" (Valleylab), 3 m length

Cleaning Instruments



20.45.4545

Cleaning pistol complete with 8 adaptors

Bipolar Coagulation Forceps U type

Ø 5 mm, WL 34 cm and 45 cm



Bipolar U type forceps 34 cm

20.60.4010	Bipolar coagulation forceps fluted bit round plug at end
20.60.4011	Bipolar coagulation forceps rippled bar with round plug at end
20.60.4012	Bipolar tubal coagulation forceps with tongue type bit with round plug at end
20.60.5010	Bipolar coagulation forceps fluted bit flat plug at end
20.60.5011	Bipolar coagulation forceps rippled bar with flat plug at end
20.60.5013	Bipolar tubal coagulation forceps with tongue type bit with flat plug at end




Bipolar U type forceps 45 cm

20.60.4013	Bipolar forceps fluted bit round plug at end
20.60.4014	Bipolar coagulation forceps rippled bar with round plug at end
20.60.4015	Bipolar tubal coagulation forceps with tongue type bit with round plug at end
20.60.5014	Bipolar coagulation forceps fluted bit flat plug at end
20.60.5015	Bipolar coagulation forceps rippled bar with flat plug at end
20.60.5016	Bipolar tubal coagulation forceps with tongue type bit with flat plug at end



Bipolar Coagulation Forceps - Powergrip

Detachable & Rotatable Ø 5 mm



20.60.4100	Handle
20.60.4101	Sheath 34 cm

This unique Bipolar instrument is used for grasping, cutting and coagulation and offers important further advantages -

- Exchangeable electrode inserts for different applications.
- Easy handling.
- Excellent cleaning and ease of repair.

Quick exchange of electrode inserts -

The different inserts can be exchanged in only a few seconds. The electrode is screwed into the sheath inserted in the handle and locked in place by means of a click-stop. At the same time the connection to the integrated electric connector is established.

Hygiene and economic efficiency -

The power grip instrument can be dismantled in three parts. This enables thorough and reliable cleaning and facilitates the execution of repair work or the replacement of defective parts.

Bipolar Inserts

For coagulation forceps - Powergrip

20.60.4110

Large grasping insert 34 cm



20.60.4111

Bipolar insert delicate curved 34 cm



20.60.4112

Bipolar grasper micro france 34 cm



20.60.4113

Bipolar grasping insert curved 34 cm



20.60.4114

Scissor insert, curved



20.60.4115

Maryland dissector, curved



20.60.4116

Claw insert



Orbitaris Bipolar Forceps



- **ORBITARIS bipolar forceps with rotating handle.**
This new innovative design provides the surgeon with comfortable ergonomic handling in all working positions. The handles rotate independently of the shaft, so the instrument can be controlled in a relaxed natural way, reducing fatigue and enabling precise handling and control during surgery.
- **Simple one handed operation**
The extended flanks on the electrode shaft, give the surgeon enhanced single finger control over the full 360° rotation available, resulting in precise control of the electrodes position. The shaft and the electrode are inserted and locked into the handle with a click stop as in the well proven POWERGRIP bipolar forceps. The electrode inserts of the POWERGRIP are fully compatible with the ORBITARIS handle resulting in a wide range of proven electrodes being available for use. This also gives POWERGRIP customers wishing to upgrade to the enhanced ergonomic handling of ORBITARIS an economic route.
- **Convenient locking device**
The ORBITARIS can be supplied with an ergonomically positioned locking device that can be operated by both right and left handed users enables the surgeon to lock the rotating handle in the chosen position. Just as easily, the locking device can be disengaged in order to move the rotating handle to a new position.
- **Rotating handle**
The handle can be rotated through 360°, so the surgeons ideal ergonomic handling position can be achieved. The surgeon and patient requirement for function and safety has been greatly enhanced, over standard handles, to include surgeon comfort during the demands of minimally invasive surgery.

20.60.0000	Orbitaris Handle
20.60.0004	Orbitaris Handle with locking device
20.60.0034	Orbitaris sheath 34 cm

Bipolar Inserts

For Orbitaris Bipolar Forceps



New inserts "P-Coat"

For these inserts a new type of insulation coating is used "P-Coat". This ultra-thin coating enables the insulation of the joints and the external surfaces of the jaws without any negative impact on their mechanical rigidity. This allows the surgeon to work very safely by ensuring that coagulation current is only applied to the chosen tissue and not the surrounding area. All P-Coat inserts are fully compatible with both the Orbitaris and the Powergrip bipolar instruments.

Orbitris & Powergrip Bipolar P-Coat inserts

20.60.8010

Large Grasping insert 34 cm with P-Coat insulation



20.60.8011

Grasping insert 34 cm fenestrated, curved with P-Coat insulation



20.60.8012

Scissor insert curved 34 cm with P-Coat insulation



20.60.8013

Dissector Maryland insert curved 34 cm with P-Coat insulation



Bipolar Cables

3 m length


	Bipolar-Units	Instruments
20.67.2003	 Erbe/Wisap/Storz	 flat plug
20.67.2103	 GIMMI / Martin / Berchtold	 flat plug
20.67.2203	 Valleylab	 flat plug

SLIMLINE Bipolar

Micro Coagulation Forceps

Reusable instruments for use in laparoscopy with a shaft Ø of 3 mm, intended to be used for fine coagulation on of tissue and small vessels using bipolar current.

The SLIMLINE Micro Coagulation Forceps can be easily detached for cleaning. The inserts can be exchanged for different applications.



SLIMLINE Bipolar Micro Coagulation Forceps

20.60.3000

3 mm, consisting of forceps body, tube, tube guide, WL 34 cm, without insert with standard handle

20.60.3001

3 mm, consisting of forceps body, tube, tube guide, WL 34 cm, without insert with ring handle

20.60.3010

Insert Rippled Bar, Jaw 3.0 mm

20.60.3011

Insert Fluted Bit, Jaw 1.5 mm



Bipolar Cables

For SLIMLINE Bipolar Micro Coagulation Forceps



20.67.3050

Cable for slimline forceps, generator side 2 banana plug, length 3 m

20.67.3055

Cable for slimline forceps, generator side Valleylab, length 3 m



20.67.3060

Cable for slimline forceps, generator side Martin, length 3 m



20.67.3065

Cable for slimline forceps, generator side Erbe, length 3 m



Bipolar Coagulation Electrodes

Ø 5 mm, 34 cm long



Bipolar Coagulation Electrodes for connection to the bipolar output of electrosurgical generators with dedicated parameters for bipolar coagulation. The Bipolar Coagulation Electrodes for Endoscopy provide the ability to dissect and coagulate tissue. The program includes instruments for coagulation of larger areas as well as of small selective areas. Due to the special design of the electrodes, they are able to withstand big mechanical stress for safe handling with consistently good coagulation results.

- A slightly conical shaft end makes insertion into the trocar sleeve easier.
- Use of bipolar technology and the optimized configuration of the electrical poles, guarantees safe and gentle work with an exactly defined area of coagulation.
- Constantly good coagulation results are achieved due to the insulation between the electrodes, always the same distance. Moreover, the electrodes are easier and more reliable to clean.
- The ergonomically designed handle allows secure holding and easy rotation during work. The standard bipolar connector plug will fit conventional bipolar cables that are readily available.

20.60.5031	Bipolar electrode, Ø 5 mm, 34 cm long, J hook
-------------------	---

20.60.5032	Bipolar electrode, Ø 5 mm, 34 cm long, L hook
-------------------	---

20.60.5033	Bipolar electrode, Ø 5 mm, 34 cm long, spatula
-------------------	--

20.60.5034	Bipolar electrode, Ø 5 mm, 34 cm long, dome
-------------------	---

20.60.5035	Bipolar electrode, Ø 5 mm, 34 cm long, plate
-------------------	--



SLIMLINE Bipolar Hooks and Needle


Ø 3 mm, WL 34 cm



Different bipolar hook-electrodes and a delicate needle electrodes provide the ability of both, a standard coagulation and a stitch coagulation. Connection with a standard bipolar forceps cable is with flat plug.


20.60.3020

SLIMLINE Bipolar electrode, L hook




20.60.3030

SLIMLINE Bipolar electrode, V hook



20.60.3040

SLIMLINE Bipolar electrode, Needle electrode



Bariatric Laparoscopy

Maxer offers following instruments for Bariatric Laparoscopy


- Laparoscopes 10 mm & 5.5 mm in longer length
- Trocars with longer tubes in 6 mm & 11 mm size
- Scissors, Graspers & Dissectors in length of 450 mm. Please add "B" to the article number while sending requests.
- The Graspers can also be supplied with metal shaft for higher bending resistance. Please add "M" to the article number while sending requests.
- Other instruments such as Veress Needles, Needle Holders Suction Irrigation Tubes, Hook Electrodes, Clip Applicators can be produced in longer working length.
- Flexible Liver Retractor along with stand to avoid damage to Liver and fatigue to Surgeon.

Single Incision Laparoscopic Surgery (SILS)

Maxer offers 5.5 mm laparoscope in 424 mm length for single port Laparoscopy. The extra length causes the camera head to shift away from instrument handles and allows more space for the operating surgeon to ease up maneuvering.



Graspers, Dissectors & Scissors



Special features of the Instruments are :

Instruments consists of the following components :

- Detachable Handle
- Detachable Tube
- Instrument insert
- Easy french ratchet

Detach ability :

- Instrument consists of 3 components – therefore easy disassembly
- Modular composition facilitates free interchange of inserts and handles
- No costly accumulation of instruments
- Optimal cleaning and sterilisation

Handling :

- Free Rotation, 360° rotation ensures more versatility
- Oblique HF connection (gold plated)
- Ergonomic, easy and precise single-hand use
- Available in Ø 3.0 mm, Ø 5.0 mm, and Ø 10.0 mm
- Light weight – therefore greater comfort

Cost effective :

- Handles with and without ratchet mechanism
- Ratchet mechanism able to be engaged
- Fully autoclavable
- Components separately available
- Only damaged component not complete instrument, needs to be replaced
- Versatility means fewer instruments are required

Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

Cobra grasping forceps

20.25.1501	20.25.1601	20.25.1701
-------------------	-------------------	-------------------

Cobra grasping forceps, double action

20.25.1502	20.25.1602	20.25.1702
-------------------	-------------------	-------------------

Allis grasping forceps

20.25.1503	20.25.1603	20.25.1703
-------------------	-------------------	-------------------

Allis grasping forceps, fenestrated

20.25.1504	20.25.1604	20.25.1704
-------------------	-------------------	-------------------

Allis grasping forceps, small

20.25.1505	20.25.1605	20.25.1705
-------------------	-------------------	-------------------

Allis grasping forceps, fenestrated, long jaw

20.25.1506	20.25.1606	20.25.1706
-------------------	-------------------	-------------------

Grasping forceps with teeth, double action

20.25.1507	20.25.1607	20.25.1707
-------------------	-------------------	-------------------

Grasping forceps, 3 : 3 teeth, double action

20.25.1508	20.25.1608	20.25.1708
-------------------	-------------------	-------------------

Grasping forceps, 2 : 1 teeth

20.25.1509	20.25.1609	20.25.1709
-------------------	-------------------	-------------------



Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

Debaquey grasping forceps, jaw length 20 mm, double action

20.25.1510	20.25.1610	20.25.1710
-------------------	-------------------	-------------------

Grasping forceps, blunt

20.25.1511	20.25.1611	20.25.1711
-------------------	-------------------	-------------------

Dissecting and grasping forceps, blunt

20.25.1512	20.25.1612	20.25.1712
-------------------	-------------------	-------------------

Dissecting and grasping forceps, with needle, blunt

20.25.1513	20.25.1613	20.25.1713
-------------------	-------------------	-------------------

Dolphin dissecting and grasping forceps, double action

20.25.1514	20.25.1614	20.25.1714
-------------------	-------------------	-------------------

Dolphin dissecting and grasping forceps, with spoon, double action

20.25.1515	20.25.1615	20.25.1715
-------------------	-------------------	-------------------

Dissecting and grasping forceps, with spoon, tapered end, double action

20.25.1516	20.25.1616	20.25.1716
-------------------	-------------------	-------------------

Grasping forceps, fenestrated, double action

20.25.1517	20.25.1617	20.25.1717
-------------------	-------------------	-------------------

Duck bill dissecting and grasping forceps

20.25.1518	20.25.1618	20.25.1718
-------------------	-------------------	-------------------



Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

Cleeve duck bill dissecting and grasping forceps, double action

20.25.1519	20.25.1619	20.25.1719
-------------------	-------------------	-------------------

Grasping forceps double action, flat

20.25.1520	20.25.1620	20.25.1720
-------------------	-------------------	-------------------

Needle nose grasping forceps

20.25.1521	20.25.1621	20.25.1721
-------------------	-------------------	-------------------

Needle nose grasping forceps, double action

20.25.1522	20.25.1622	20.25.1722
-------------------	-------------------	-------------------

Maryland dissecting forceps, slightly curved, double action

20.25.1523	20.25.1623	20.25.1723
-------------------	-------------------	-------------------

Kelly grasping forceps, double action

20.25.1524	20.25.1624	20.25.1724
-------------------	-------------------	-------------------

Grasping forceps, fenestrated, double action, delicate model

20.25.1525	20.25.1625	20.25.1725
-------------------	-------------------	-------------------

Fabzang, small, double action

20.25.1526	20.25.1626	20.25.1726
-------------------	-------------------	-------------------

Grasping forceps, fenestrated, strong, double action

20.25.1527	20.25.1627	20.25.1727
-------------------	-------------------	-------------------



Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

Grasping forceps, fenestrated, strong, double action

20.25.1527	20.25.1627	20.25.1727
-------------------	-------------------	-------------------

Mixer dissecting forceps, 90° angled, double action

20.25.1528	20.25.1628	20.25.1728
-------------------	-------------------	-------------------

Grasping forceps, 1 : 2 teeth, double action

20.25.1529	20.25.1629	20.25.1729
-------------------	-------------------	-------------------

Dissecting forceps, not serrated, double action

20.25.1530	20.25.1630	20.25.1730
-------------------	-------------------	-------------------

Grasping forceps, double action

20.25.1531	20.25.1631	20.25.1731
-------------------	-------------------	-------------------

Grasping forceps, fenestrated, strong, double action

-	20.25.1632	20.25.1732
---	-------------------	-------------------

Grasping forceps, atraumatic, double action

20.25.1533	20.25.1633	20.25.1733
-------------------	-------------------	-------------------

Grasping forceps, double action

20.25.1534	20.25.1634	20.25.1734
-------------------	-------------------	-------------------

Bowel forceps, double action

20.25.1535	20.25.1635	20.25.1735
-------------------	-------------------	-------------------



Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

Tube grasping forceps, double action

20.25.1536	20.25.1636	20.25.1736
-------------------	-------------------	-------------------

Dissecting and grasping forceps, blunt, with spoon, double action

20.25.1537	20.25.1637	20.25.1737
-------------------	-------------------	-------------------

Heavy duty (fundus) grasping forceps, double action

20.25.1538	20.25.1638	20.25.1738
-------------------	-------------------	-------------------

Alligator grasping forceps, double action

20.25.1539	20.25.1639	20.25.1739
-------------------	-------------------	-------------------

Babcock grasping forceps, cross serrated

20.25.1540	20.25.1640	20.25.1740
-------------------	-------------------	-------------------

Babcock grasping forceps, jaw length 30 mm, double action

20.25.1541	20.25.1641	20.25.1741
-------------------	-------------------	-------------------

Dissecting forceps, double action, flat

20.25.1542	20.25.1642	20.25.1742
-------------------	-------------------	-------------------

Dissecting forceps, delicate pattern, flat

20.25.1543	20.25.1643	20.25.1743
-------------------	-------------------	-------------------

Schertl grasping forceps, fenestrated, double action

20.25.1544	20.25.1644	20.25.1744
-------------------	-------------------	-------------------



Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

Michel atraumatic dissecting forceps

20.25.1545	20.25.1645	20.25.1745
-------------------	-------------------	-------------------

Michel atraumatic dissecting forceps, double action

20.25.1546	20.25.1646	20.25.1746
-------------------	-------------------	-------------------

Grasping forceps with tungsten carbide inserts, double action

20.25.1547	20.25.1647	20.25.1747
-------------------	-------------------	-------------------

Micro grasping forceps

20.25.1548	20.25.1648	20.25.1748
-------------------	-------------------	-------------------

Dorsey intestinal forceps, jaw length 40 mm, double action

20.25.1549	20.25.1649	20.25.1749
-------------------	-------------------	-------------------

Claw grasping forceps

20.25.1550	20.25.1650	20.25.1750
-------------------	-------------------	-------------------

Allis grasping forceps, jaw length 30 mm, double action

20.25.1551	20.25.1651	20.25.1751
-------------------	-------------------	-------------------

Hunter bowel grasping forceps, double action

20.25.1552	20.25.1652	20.25.1752
-------------------	-------------------	-------------------



Design changes may take place without prior notice.

L - 63

Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

De bailey grasping forceps, jaw length 30 mm, double action

20.25.1553	20.25.1653	20.25.1753
-------------------	-------------------	-------------------

Biopsy and grasping forceps, double action

20.25.1554	20.25.1654	20.25.1754
-------------------	-------------------	-------------------

Biopsy forceps, spoon jaw with needle, double action

20.25.1555	20.25.1655	20.25.1755
-------------------	-------------------	-------------------

Biopsy forceps, spoon jaw, double action

20.25.1556	20.25.1656	20.25.1756
-------------------	-------------------	-------------------

Biopsy forceps, spoon jaw with needle

20.25.1557	20.25.1657	20.25.1757
-------------------	-------------------	-------------------

Biopsy forceps, spoon jaw

20.25.1558	20.25.1658	20.25.1758
-------------------	-------------------	-------------------

Blakesley biopsy forceps

20.25.1559	20.25.1659	20.25.1759
-------------------	-------------------	-------------------

Biopsy forceps, cut - thru type

20.25.1560	20.25.1660	20.25.1760
-------------------	-------------------	-------------------



Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

Frangenheim biopsy forceps, cut - thru type

20.25.1561	20.25.1661	20.25.1761
-------------------	-------------------	-------------------

Mini biopsy forceps, spoon - jaw

20.25.1562	20.25.1662	20.25.1762
-------------------	-------------------	-------------------

Peritoneal Scissors, toothed, blunt, straight

20.25.5501	20.25.5601	-
-------------------	-------------------	---

Peritoneal scissors, toothed, sharp tips, straight

20.25.5502	20.25.5602	-
-------------------	-------------------	---

Hook scissors

20.25.5503	20.25.5603	-
-------------------	-------------------	---

Mini Hook scissors

20.25.5504	20.25.5604	-
-------------------	-------------------	---

Metzenbaum Scissors, straight, double action

20.25.5505	20.25.5605	-
-------------------	-------------------	---

Metzenbaum EndoScissors, strong curved, double action

20.25.5506	20.25.5606	-
-------------------	-------------------	---

Metzenbaum Scissors, slightly curved, double action

20.25.5507	20.25.5607	-
-------------------	-------------------	---



Graspers, Dissectors & Scissors

Ø 5 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010

Metzenbaum Scissors, straight, long blade, double action

20.25.5508	20.25.5608	-
-------------------	-------------------	---

Metzenbaum Scissors, curved, long blade, double action

20.25.5509	20.25.5609	-
-------------------	-------------------	---

Metzenbaum Scissors, curved, toothed, long blade, double action

20.25.5510	20.25.5610	-
-------------------	-------------------	---

Baby Metzenbaum Scissors, curved, double action

20.25.5511	20.25.5611	-
-------------------	-------------------	---

Metzenbaum Scissors, straight, short, pointed, double action

20.25.5512	20.25.5612	-
-------------------	-------------------	---

Mini-Metzenbaum Scissors, curved, double action

20.25.5513	20.25.5613	-
-------------------	-------------------	---

Mini-Metzenbaum Scissors, toothed, curved, double action

20.25.5514	20.25.5614	-
-------------------	-------------------	---

Micro scissors, straight, pointed

20.25.5515	20.25.5615	-
-------------------	-------------------	---

Micro scissors, curved, pointed

20.25.5516	20.25.5616	-
-------------------	-------------------	---



Graspers, Dissectors & Scissors

Ø 10 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 10 mm with knob
	20.25.0001	20.25.0002	20.25.0310

Claw Grasping forceps

20.30.0501	20.30.0601	20.30.0701
-------------------	-------------------	-------------------

Claw Grasping forceps, double action

20.30.0502	20.30.0602	20.30.0702
-------------------	-------------------	-------------------

OP Scissors

20.30.0503	20.30.0603	-
-------------------	-------------------	---

OP Spoon forceps

20.30.0504	20.30.0604	20.30.0704
-------------------	-------------------	-------------------

Allis Grasping forceps, double action

20.30.0505	20.30.0605	20.30.0705
-------------------	-------------------	-------------------

Duval Grasping forceps, double action

20.30.0506	20.30.0606	20.30.0706
-------------------	-------------------	-------------------

Babcock Grasping forceps, double action

20.30.0507	20.30.0607	20.30.0707
-------------------	-------------------	-------------------

Marc grasping forceps, double action

20.30.0508	20.30.0608	20.30.0708
-------------------	-------------------	-------------------



Graspers, Dissectors & Scissors

Ø 10 mm, WL 33 cm, 360° Rotatable and detachable

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 10 mm with knob
	20.25.0001	20.25.0002	20.25.0310

Pennington grasping forceps, double action

20.30.0509	20.30.0609	20.30.0709
-------------------	-------------------	-------------------

Debaquey intestinal forceps, jaw length 40 mm, double action

20.30.0510	20.30.0610	20.30.0710
-------------------	-------------------	-------------------

Koeln grasping forceps, double action

20.30.0511	20.30.0611	20.30.0711
-------------------	-------------------	-------------------

Alligator grasping forceps, double action

20.30.0512	20.30.0612	20.30.0712
-------------------	-------------------	-------------------

Paddle babcock grasping forceps, double action

20.30.0513	20.30.0613	20.30.0713
-------------------	-------------------	-------------------

Mixer dissecting forceps, 90° angled, double action

20.30.0514	20.30.0614	20.30.0714
-------------------	-------------------	-------------------

Maryland dissecting forceps, slightly curved, double action

20.30.0515	20.30.0615	20.30.0715
-------------------	-------------------	-------------------

Maryland dissecting forceps, strongly curved, double action

20.30.0516	20.30.0616	20.30.0716
-------------------	-------------------	-------------------



Pediatric Laparoscopy

Graspers, Dissectors & Scissors - 360° Rotatable and detachable

Pediatric Ø 3 mm, WL 20 cm, Instruments P at end for inserts and complete instruments article no.

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 3 mm with knob
	20.25.0001	20.25.0002	20.25.0010P

Example : Cobra Grasping forceps : Pediatric Ø 3 mm, WL 20 cm

20.25.1501P + 20.25.0001 + 20.25.0010P = 20.25.1601P

20.25.1501P + 20.25.0002 + 20.25.0010P = 20.25.1701P

Pediatric Ø 5 mm, WL 25 cm, Instruments S at end for inserts and complete instruments article no.

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010S

Example : Cobra Grasping forceps : Pediatric Ø 5 mm, WL 25 cm

20.25.1501S + 20.25.0001 + 20.25.0010S = 20.25.1601S

20.25.1501S + 20.25.0002 + 20.25.0010S = 20.25.1701S

Bariatric Laparoscopy

Graspers, Dissectors & Scissors - 360° Rotatable and detachable

Bariatric Ø 5 mm, WL 45 cm, Instruments B at end for inserts and complete instruments article no.

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 5 mm with knob
	20.25.0001	20.25.0002	20.25.0010B

Example : Cobra Grasping forceps : Bariatric Ø 5 mm, WL 45 cm

20.25.1501B + 20.25.0001 + 20.25.0010B = 20.25.1601B

20.25.1501B + 20.25.0002 + 20.25.0010B = 20.25.1701B

Bariatric Ø 10 mm, WL 45 cm, Instruments B at end for inserts and complete instruments article no.

			
Insert only	Handle without Ratchet Plastic handle with insulated shaft, without ratchet, with HF connection, rotatable 360°	Handle with Ratchet Plastic handle with insulated shaft, with ratchet, with HF connection, rotatable 360°	Insulating detachable tube Ø 10 mm with knob
	20.25.0001	20.25.0002	20.25.0310B

Example : Claw Grasping forceps : Bariatric Ø 10 mm, WL 45 cm

20.30.0501B + 20.25.0001 + 20.25.0310B = 20.30.0601B

20.30.0501B + 20.25.0002 + 20.25.0310B = 20.30.0701B

GERMANY
LAPAROSCOPY[illegible]

GERMANY
LAPAROSCOPY[illegible]

GERMANY
LAPAROSCOPY[illegible]



GERMANY

UNITS

LAPAROSCOPY

UROLOGY

HYSTEROSCOPY

ARTHROSCOPY

ENT

PEDIATRIC

SPINE ENDOSCOPY

Distributed by:

Maxer Medizintechnik GmbH

Untere Hauptstr. 34/1

78573 Wurmlingen

Germany

Registered Office: Vogesenstr.17

78549 Spaichingen

Germany

Email: info@maxerendoscopy.com

Website: www.maxerendoscopy.com