

**Maxer**  
GERMANY

ACCESSORIES



**M** ACCESSORIES  
**SYSTEM**

[maxerendoscopy.com](http://maxerendoscopy.com)





## ISO Certificate

Quality management system as per  
ISO 13485 :2012 Medical Device Certification  
for development & Production & Endoscopic  
instrument & Equipments

## CE Certificate

As per council directive 93/42/EEC  
Concerning medical devices





## About Maxer

Maxer Medizintechnik GmbH has established itself as a symbol of creativity and latest technology and is an internationally qualified supplier of world class medical instruments and devices.

Specialized in the field of Minimal Invasive Surgery, Maxer manufactures wide range of equipments such as Laparoscopy, Hysteroscopy, Urology, Pediatric Surgery, Arthroscopy, Spine Endoscopy, ENT, etc.

Maxer has been one of the foremost companies to bring out HD technology in MIS field. In 2004, Maxer brought out first endoscopy camera with Matrix Metering. In 2008, Maxer brought out HD Telescopes. In 2009, Maxer launched first Networked Endoscopy Camera System. In 2010, Maxer introduced pathbreaking 16 CH Resectoscope for Office Hysteroscopy. In 2011, Maxer launched LED Headlight with 200,000 Lux illumination, again a state of art in this segment. In 2012, Maxer launched LED Light sources to keep Maxer customers always at the cutting edge of technology. Maxer has also launched HD cameras with integrated Full HD recording. 3 CMOS HD cameras with 6 megapixel camera head will be launched soon.

Owners of Maxer have been working in Medical market for more than 20 years and have installed thousands of sets with high level of customer satisfaction. This background has given Maxer profound credibility in the field of endoscopy.

Due to vibrant co-operation with end users, Maxer is able to continuously improve upon the product range and also able to meet special modification demands. Maxer believes that growth and innovation go hand-in-hand.




Maxer sells complete system with a focus on product quality and customers satisfaction. Maxer provides sales and service training to it's distributor in order to be successful.

Maxer believes that sound ethics with professionalism will only give success to the company. Maxer's growth and success are based on long term partnerships with customers, distributors. Maxer currently exports it's products to Europe, Asia, Africa.

Maxer is headed by Mr. Giovanni Miranda in the capacity as a Managing Director.

How to reach to Maxer - Maxer is located in Würlingen (sales office) within the medical instruments manufacturing region around Tuttlingen in South Germany. Maxer also has registered office at Spaichingen.







## INDEX

Basic Trolley	T - 4 To 5
Easy Trolley	T - 6
Endo Trolley	T - 7
Endo Arm Trolley	T - 8 To 9
Maxer Trolley - Accessories	T - 10 To 11
Stainless Steel Max Cart	T - 12
Sterilization Wire Baskets	T - 13
Sterilization Trays	T - 14 To 16
Sterilization Containers & Endoscope Holder	T - 17 To 18
Cleaning Instruments	T - 19



## Basic Trolley



**10.82.3001**

Basic trolley with 3 shelves, drawer, handle, Socket strip 6 way, 4 twin swivel wheels lockable, total height 1117 mm, width 500 mm

**10.82.3001W**

Basic trolley with 3 shelves, drawer, handle, Socket strip 6 way, 4 twin swivel wheels lockable, total height 1117 mm, width 600 mm



**10.82.3002**

Basic trolley with 4 shelves, drawer, handle, Socket strip 6 way, 4 twin swivel wheels lockable, total height 1117 mm, width 500 mm

**10.82.3002W**

Basic trolley with 4 shelves, drawer, handle, Socket strip 6 way, 4 twin swivel wheels lockable, total height 1117 mm, width 600 mm



## Basic Trolley



<b>10.82.3003</b>	Basic trolley with 5 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, total height 1437 mm, width 500 mm
<b>10.82.3003W</b>	Basic trolley with 5 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, total height 1437 mm, width 600 mm



## Easy Trolley



**10.82.3011**

Easy trolley with 3 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, Monitor holder on cross piece, total height 1117 mm, width 500 mm

**10.82.3011W**

Easy trolley with 3 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, Monitor holder on cross piece, total height 1117 mm, width 600 mm



**10.82.3012**

Easy trolley with 4 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, Monitor holder on cross piece, total height 1437 mm, width 500 mm

**10.82.3012W**

Easy trolley with 4 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, Monitor holder on cross piece, total height 1437 mm, width 600 mm



## Endo Trolley



**10.82.3021**

Endo trolley with 3 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, Monitor holder on Shelf, total height 1117 mm, width 500 mm

**10.82.3021W**

Endo trolley with 3 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, Monitor holder on Shelf, total height 1117 mm, width 600 mm



**10.82.3022**

Endo trolley with 4 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, Monitor holder on Shelf, total height 1437 mm, width 500 mm

**10.82.3022W**

Endo trolley with 4 shelves, drawer, handle, Socket strip 6-way, 4 twin swivel wheels lockable, Monitor holder on Shelf, total height 1437 mm, width 600 mm



## Endo Arm Trolley




**10.82.3032**

Easy arm trolley with 4 shelves, drawer, handle, Socket strip 6-way, 4 swivel lockable wheels, **Monitor arm on cross piece**, with counter weights, total height 1117 mm, width 500 mm

**10.82.3032W**

Easy arm trolley with 4 shelves, drawer, handle, Socket strip 6-way, 4 swivel lockable wheels, **Monitor arm on cross piece**, with counter weights, total height 1117 mm, width 600 mm



**10.82.3033**


Easy arm trolley with 5 shelves, drawer, handle, Socket strip 6-way, 4 swivel lockable wheels, **Monitor arm on cross piece**, with counter weights, total height 1437 mm, width 500 mm

**10.82.3033W**

Easy arm trolley with 5 shelves, drawer, handle, Socket strip 6-way, 4 swivel lockable wheels, **Monitor arm on cross piece**, with counter weights, total height 1437 mm, width 600 mm



## Endo Arm Trolley




**10.82.3042**

Endo arm trolley with 4 shelves, drawer, handle, Socket strip 6-way, 4 swivel lockable wheels, **Monitor arm & monitor holder on Shelf**, with counter weights, total height 1117 mm, width 500 mm

**10.82.3042W**

Endo arm trolley with 4 shelves, drawer, handle, Socket strip 6-way, 4 swivel lockable wheels, **Monitor arm & monitor holder on Shelf**, with counter weights, total height 1117 mm, width 600 mm



**10.82.3043**

Endo arm trolley with 5 shelves, drawer, handle, Socket strip 6-way, 4 swivel lockable wheels, **Monitor arm & monitor holder on Shelf**, with counter weights, total height 1437 mm, width 500 mm

**10.82.3043W**

Endo arm trolley with 5 shelves, drawer, handle, Socket strip 6-way, 4 swivel lockable wheels, **Monitor arm & monitor holder on Shelf**, with counter weights, total height 1437 mm, width 600 mm

Design changes may take place without prior notice.

T - 9



**Maxer Trolley - Accessories**

<b>10.82.1350</b>	Isolation transformer 660 VA
<b>10.82.1351</b>	Isolation transformer with display 660 VA
<b>10.82.1352</b>	Isolation transformer 1000 VA
<b>10.82.1353</b>	Isolation transformer with display 1000 VA
<b>10.82.1354</b>	Isolation transformer 1600 VA
<b>10.82.1355</b>	Isolation transformer with display 1600 VA



## Maxer Trolley - Accessories



**10.82.1004**

Camera holder with clamp



**10.82.1007**

Hose holder



**10.82.1005**

Gas bottle holder  
(Capacity 20 kg)



**10.82.1008**

Cable holder



**10.82.1006**

Foot switch cradle



**10.82.1009**

Standard runner 10 x 25 mm,  
L390



**Maxer Trolley - Max Cart**




**10.80.3000**

Trolley complete in Stainless steel with 4 shelves, one drawer & cable socket up to 6 cables



## Sterilization Accessories

### Wire basket




**15.10.0021**

Wire basket stainless steel with cover for autoclaving 310 x 80 x 60 mm with silicone bands for 185 mm long Optics

**15.10.0011**

Wire basket stainless steel with cover for autoclaving 450 x 80 x 60 mm with silicone bands for 330 mm long Optics



**15.10.0075**

Wire basket stainless steel for autoclaving 640 x 150 x 77 mm with silicon holders and lid for URS WL 425 mm


**15.10.0055**

Wire basket stainless steel for autoclaving 670 x 80 x 55 mm with silicon holders and lid for Bronchoscopes WL 490 mm




## Sterilization Accessories

### Stainless Steel Trays



**15.13.0010**

Stainless Steel Tray for autoclaving 484 x 252 x 65 mm (without instruments)




**15.13.0012**

Stainless Steel Tray for autoclaving 484 x 252 x 65 mm (without instruments)




## Sterilization Accessories

### Stainless Steel Trays



**15.13.0011**

Stainless Steel Tray for autoclaving 484 x 252 x 65 mm (without instruments)



**15.13.0013**

Stainless Steel Tray for autoclaving 539 x 252 x 95 mm (without instruments)



## Sterilization Accessories

### Stainless Steel Trays




**60.50.1010**







Spine Instruments sterilization tray (without instruments)









## Sterilization Accessories

Aluminium container




	Aluminium container 580x280x200 mm for the wire baskets	
<b>15.15.1000</b>	Grey	
<b>15.15.1000YL</b>	Yellow	
<b>15.15.1000RD</b>	Red	
<b>15.15.1000BL</b>	Blue	
<b>15.15.1000GN</b>	Green	
<b>15.15.1000BK</b>	Black	

	Aluminium container 580x280x260 mm for the wire baskets	
<b>15.15.1001</b>	Grey	
<b>15.15.1001YL</b>	Yellow	
<b>15.15.1001RD</b>	Red	
<b>15.15.1001BL</b>	Blue	
<b>15.15.1001GN</b>	Green	
<b>15.15.1001BK</b>	Black	



**Sterilization Accessories & Endoscope Holder**



**15.15.1100**

Endoscopy instrument stand 485 x 235 x 160 mm with silicon holder bar for 10 instruments

**Endoscope Holder**



**12.03.0010**


Endoscope and liver retractor holder

**12.03.0011**

Clamp for holder



**Cleaning Instruments**



**20.45.4545**

Cleaning pistol complete with 8 adaptors



[illegible]





UNITS

LAPAROSCOPY

UROLOGY

HYSTEROSCOPY

ARTHROSCOPY

ENT

PEDIATRIC

SPINE ENDOSCOPY

Distributed by:

**Maxer Medizintechnik GmbH**

Untere Hauptstr. 34/1

78573 Wurmlingen

Germany

Registered Office: Vogesenstr.17

78549 Spaichingen

Germany

Email: [info@maxerendoscopy.com](mailto:info@maxerendoscopy.com)

Website: [www.maxerendoscopy.com](http://www.maxerendoscopy.com)

VERSION 15H28